


DISPOSABLES AND ACCESSORIES CATALOGUE



www.laborie.com

URODYNAMICS - BIOFEEDBACK / STIMULATION - ULTRASOUND - ANORECTAL MANOMETRY

LABORIE MEDICAL TECHNOLOGIES

Quality, Innovation, Service and Satisfaction

Dear Valued Customer,

We are pleased to introduce our new Disposables & Accessories Catalogue which features information about our complete range of disposable and reusable accessories.

A concise list of all products can be found at the end for quick reference; all of our products and their descriptions can be found under specific applications as listed below.

- Urodynamics
- Electromyography
- Biofeedback / Stimulation
- Anorectal Manometry
- Ultrasound

In 2007 we celebrated our 40th anniversary in the field of Urodynamics; we are privileged to partner with numerous leading physicians and specialists in the field. This catalogue outlines the safe and cost-effective products which were developed in response to their requirements.

Our company continues to strive for quality, innovation, service, and satisfaction. From ergonomic designs to accurate package labelling, our products ensure confidence from the time a package is opened until the end of an examination. Quality and satisfaction are among our trademarks and can be seen throughout all Laborie products.

It is our hope that you will find this catalogue useful; please do not hesitate to contact your local representative for detailed information.

Sincerely Yours,

Ray Laborie

Chairman



TABLE OF CONTENTS


■ URODYNAMICS	
Catheters / Fluid Filled	1-2
Accessories / Fluid Filled	3
Catheters / Air Charged	4
Accessories / Air Charged	5
Tubing	6
■ ELECTROMYOGRAPHY	
Patches and Leads	7
■ BIOFEEDBACK / STIMULATION	
Probes	8
Patch Electrodes	9
Manometry Probes and Accessories	10
■ ANORECTAL MANOMETRY	
Catheters and Accessories	11
■ ULTRASOUND	
Covers and Gel	12
■ MISCELLANEOUS	
Funnels, Beakers and Commodes	13
■ INDEX	
Reference Guide	14-15
■ HARDWARE	
Systems	16-18




Some class II devices advertised may not be
licensed for sale in accordance with Canadian law.
Please check with your local representative.

URODYNAMICS

Catheters / Fluid Filled



Twin Lumen



Triple Lumen

For measuring bladder and urethral pressures, Laborie's urodynamic catheters made from a custom Pebax blend, offer a variety of features.


- Optimized flexibility ensures patient comfort during catheterization
- Luer lock connectors with ergonomic wings for easy connection
- Universal color coded channels for quick and easy identification
- Lumens designed to reduce risk of channel crosstalking (filling lumen included)
- The pressure hole is cut with advanced technology to maximize patient comfort

CATHETERS


Item no.	Description	Shaft Size	Length	Quantity / Box
CAT206	Twin Lumen	6 Fr	40cm (15.7in)	10
CAT208	Twin Lumen	8 Fr	40cm (15.7in)	10
CAT307	Triple Lumen*	7 Fr	40cm (15.7in)	10

* With graduated mark at 1cm, from 9cm, to 25cm, and a radiopaque mark at 6cm, from the distal tip.

CAT206 / CAT208 Twin Lumen




CAT307 Triple Lumen



URODYNAMICS Catheters / Fluid Filled

For abdominal pressure monitoring during Cystometry (CMG) studies, these catheters are optimally designed.

- Convenient “through and through” flushing action for ease of preparation
- Requires only 1-1.5mL of filling ensuring patient comfort
- Wing-tipped luer lock connectors ensure ease of use and a secure connection
- Latex and Silicone options available
- Pediatric option available




CATHETERS INTERNATIONAL

Item no.	Description	Shaft Size	Balloon Size	Length	Quantity / Box
CAT411	Latex Rectal Balloon (Adult)	10 Fr	5cc	40cm (15.7in)	10
CAT506	Silicone Rectal Balloon (Pediatric)	8 Fr	2.5cc	40cm (15.7in)	10
CAT509	Silicone Rectal Balloon (Adult)	10 Fr	5cc	40cm (15.7in)	10

CATHETERS USA

Item no.	Description	Shaft Size	Balloon Size	Length	Quantity / Box
CAT410	Latex Rectal Balloon (Adult)	10 Fr	5cc	40cm (15.7in)	10
CAT508	Silicone Rectal Balloon (Pediatric)	8 Fr	2.5cc	40cm (15.7in)	10
CAT510	Silicone Rectal Balloon (Adult)	10 Fr	5cc	40cm (15.7in)	10
CAT108	Silicone Rectal Balloon (Adult)	14 Fr	5cc	38cm (14.9in)	10

Rectal Balloon Catheter



URODYNAMICS

Accessories / Fluid Filled

ACCESSORIES

Item no.	Description	Quantity / Box
DIS102	Luer Lock Plugs	25
DIS103	Nova Dome	25 / 50
DIS203	Nova Dome with Luer Lock Plug (USA only)	25 / 50
DIS130	Transducer Cartridge with Luer Lock	25

MISCELLANEOUS

Item no.	Description	Quantity / Box
DIS105	Catheter Adapter (USA only)	10
CAT103	1 - Way Catheter Adapter (USA only)	1
DIS302	Infusion Cuff	1
DIS431	Stopcock - 3 - way Universal	50
DIS173	Graduated Urine Beaker	10



Nova Dome

DIS103

Disposable cover for accurate pressure transmission to reusable transducers. (For use with DIS102+ TUB101)



Luer Lock Plugs

DIS102

Ensures a securely closed dome for pressure containment. (For use with DIS130 + DIS103)



Pressure Cuff

DIS302

For use with fluid filled catheters during UPP Studies. This cuff provides consistent pressure throughout testing for accurate measurement of urethral pressure. Designed for use with 500mL or 1,000mL bags of fluid.




Transducer Cartridge with Luer Lock Plug

DIS130

Transducer cover for accurate pressure transmission to reusable transducers, includes luer lock plug.

URODYNAMICS

Air Charged Catheters



Air Charged Catheters offer a reliable circumferential pressure monitoring system.


- Easy to use.
- Single or Dual Sensor options
- Cost - Competitive
- Radio opaque markings available

Laborie's Global leading UDS software package offers an array of urodynamics test capabilities depending on your clinical needs. Contact your local representative for details.

CATHETERS

Item no.	Description	Length	Shaft Size	Sensor Distance	Quantity / Box
CAT875	Abdominal Catheter without Lumen Rectal / Vaginal	50cm (19.7in)	7 Fr	-	10
CAT878	Single Sensor with Lumen + Coude Tip Bladder	50cm (19.7in)	7 Fr	-	10
CAT880	Dual Sensor Catheter with Lumen Bladder / Urethral	50cm (19.7in)	7 Fr	6cm	10
CAT885	Dual Sensor Catheter with Lumen Bladder / Urethral (Radio - opaque)	50cm (19.7in)	7 Fr	6cm	10
CAT895	Single Sensor with Lumen Bladder / Urethral	50cm (19.7in)	7 Fr	-	10
CAT905	Single Sensor with Lumen Bladder / Urethral (Radio - opaque)	50cm (19.7in)	7 Fr	-	10

Air Charged Catheter



URODYNAMICS

Accessories / Air Charged



Stylet 10

DIS500

This stylet makes the catheter less flexible to aid with difficult catheter placements by insertion into the filling port of any single or dual sensor with lumen catheter. The stylet also offers the ability to bend, simulating a coude tip catheter. NOTE: after catheter is in place, this stylet must be removed to allow fluid infusion.

ACCESSORIES

Item no.	Description	Quantity / Box
DIS450	Air Charged Calibrator Chamber	1
DIS500	Stylet 10	1
DIS950	Velcro Leg Strap Stabilizer	1



Air Charged Calibrator Chamber

DIS450


Recommended for use when calibration is necessary. Provides a great range of pressure for accurate and stable calibration.



Velcro Leg Strap Stabilizer

DIS950

This strap allows the transducer cables to be secured to the patient's leg to increase ergonomics and reduce risk of catheter movement. Each transducer cable is fixed with Velcro and can be secured at any location along the length of the strap.



Pump Tubing Infusion Line

TUB500

Includes all components needed to accurately infuse fluid through urodynamics catheters.

TUBING

Item no.	Description	Length	Quantity / Box
TUB101	Measurement Tubing	150cm (59.1in)	25 / 50
TUB160	Clear Monitoring Line (USA only)	183cm (72.0in)	25
TUB219	Perfusion Line with 2 connections (International only)	50cm (19.7in)	5
TUB500	Pump Tubing Infusion Line	400cm (157.1in)	10 / 25
TUB501	Pump Head Tubing	22cm (8.7in)	5
TUB576	Extension Tubing without Stopcock	185cm (72.8in)	50 / 100



Measurement Tubing

TUB101

Pressure measurement tubing with attached stopcock.



Extension Tubing without Stopcock

TUB576

Extension tubing with male-male connectors



Pump Head Tubing

TUB501

Pump head tubing; replacement for the compressible portion only of Laborie pump tubing

ELECTROMYOGRAPHY

Patches and Leads



EMG Patch Electrodes, without Lead Cables, set of 3

ELE425

Surface electrodes measure pelvic floor / sphincter activity. Supplied in strips of three patches (two “active” and one ground) patches, they snap securely onto lead wires. Gel adhesive for secure placement. May also be used for EMG monitoring during biofeedback sessions using Urostym. For use with ELE350, ELE365, ELE370 re-useable lead wires.



EMG Patch Electrodes, with Lead Cables, set of 3

ELE428

EMG patch electrodes, with lead cables and surface electrodes to measure pelvic floor/sphincter activity. These patch electrodes have lead wires already attached, for easy, complete disposal after each patient use. Soft cloth backing and strong adhesive allow for secure placement. May also be used for EMG monitoring during biofeedback sessions using Urostym. For use with EMG cable CAB154 only.

PATCH ELECTRODES

Item no.	Description	Quantity / Box
ELE335	EMG Pediatric Gel Patch Electrodes Set of 3 (USA only)	10
ELE425	EMG Patch Electrodes, without Lead Cables, set of 3	10
ELE428	EMG Patch Electrodes, with Lead Cables, set of 3	10

LEADS


Item no.	Description	Quantity / Box
ELE200	EMG Lead Set / Shielded (Black/Red/White)	1
ELE365	EMG Lead Wire Non - Shielded 1 Green (USA only)	1
ELE370	EMG Lead Wires Non - Shielded 2 Red (USA only)	1
ELE350	EMG Lead Set, Non Shielded (2 Red, 1 Green)	1



EMG Lead Set / Shielded

ELE200

For use when EMG shielding is important. These leads insert securely into headstage. Requires CAB155 to connect to hardware.



EMG Lead Set, Non Shielded

ELE350

2 red and 1 green lead per set. For use with non-shielded headstages, these wires snap onto the EMG surface electrodes. For use with ELE425, ELE335 patch electrodes. Requires CAB154 to connect to hardware.

BIOFEEDBACK / STIMULATION Probes

VAGINAL PROBES

Item no.	Description	Length	Quantity / Box
PRB450	Vaginal Classic Probe	10cm(4in)	1
PRB085	Vaginal Atrophic Probe	8cm(3in)	1
PRB425	Vaginal Tampon Probe	5cm(2in)	1

ANAL PROBES

Item no.	Description	Length	Quantity / Box
PRB055	Anal Probe	8.3cm(3.25in)	1

All probes can be used for biofeedback and electrical stimulation



Vaginal Classic Probe

PRB450

Standard vaginal probe will deliver up to 70mA of electrical stimulation. All probes are for single patient use.



Vaginal Tampon Probe

PRB425

Vaginal tampon probe is capable of delivering up to 70mA of electrical stimulation. All probes are for single patient use.



Vaginal Atrophic Probe

PRB085

More slender vaginal probe, is capable of up to 70mA electrical stimulation. For use with vaginal atrophy. All probes are for single patient use.



Anal Probe

PRB055

For both biofeedback and electrical Stimulation. The anal probe is set to deliver up to 25mA of electrical stimulation. All probes are for single patient use.

BIOFEEDBACK / STIMULATION

Patch Electrodes

ELECTRODES

Item no.	Description	Quantity / Box
ELE428	EMG Patch Electrodes, with Lead Cables, set of 3	10
ELE425	EMG Patch Electrodes, without Lead Cables set of 3	10
ELE625	Surface Electrode / Stimulation Reusable set of 4	15



EMG Patch Electrodes, with Lead Cables, set of 3

ELE428

EMG patch electrodes, with lead cables and surface electrodes to measure pelvic floor/sphincter activity. These patch electrodes have lead wires already attached, for easy, complete disposal after each patient use. Soft cloth backing and strong adhesive allow for secure placement. May also be used for EMG monitoring during biofeedback sessions using Urostym. For use with EMG cable CAB154 only.



EMG Patch Electrodes, without Lead Cables, set of 3

ELE425

Surface electrodes measure pelvic floor / sphincter activity. Supplied in strips of three patches (two "active" and one ground) patches, they snap securely onto lead wires. Gel adhesive for secure placement. May also be used for EMG monitoring during biofeedback sessions using Urostym. For use with ELE350, ELE365, ELE370 re-useable lead wires.



Surface Electrode / Stimulation Reusable set of 4

ELE625

Dual purpose patches will monitor EMG activity. Useful on pelvic and abdominal muscles. (Reusable on same patient by storing in enclosed plastic film for next use)

BIOFEEDBACK / STIMULATION

Manometry Probes and Accessories



Graduated Urine Beaker

DIS173

For use during Uroflow or Pressure Flow Studies, this beaker has convenient mL markers for confirmation of fluid volumes. Designed to fit perfectly on Laborie Uroflowmeters.



Nova Dome

DIS103

Disposable cover for accurate pressure transmission to reusable transducers. (For use with DIS102, TUB101 + CAB908)

ACCESSORIES

Item no.	Description	Quantity / Box
DIS173	Graduated Urine Beaker	10
DIS203	Nova Dome with Luer Lock Plug (USA only)	25 / 50
DIS102	Luer Lock Plugs	25
DIS103	Nova Dome	25 / 50

PROBES AND TUBING

Item no.	Description	Length	Quantity / Box
PRB113	Vaginal Manometry Probe	11cm(5.25in)	1
PRB080	Vaginal Pressure Tubing	77cm(30in)	1
PRB108	Anal Manometry Probe	6cm(2.25in)	1
PRB078	Anal Pressure Tubing	122cm(48in)	1



Vaginal Manometry Probe

PRB113

Medical grade silicone elastomer (high biocompatibility). Latex-free, micro-porous, supple, air filled probe. Air is inserted / removed with a syringe to measure pressure. (PRB080) only.



Vaginal Pressure Tubing

PRB080

For use with Vaginal Manometry Probe (PRB113) only.



Anal Manometry Probe

PRB108

Medical grade silicone elastomer (high biocompatibility). Latex-free, micro-porous, supple, air filled probe. (PRB078) only.



Anal Pressure Tubing

PRB078

Pressure tubing for use with Anal Manometry probe (PRB108) only.

ANORECTAL MANOMETRY

Catheters and Accessories

Laborie's Anorectal Manometry System is capable of functioning with a variety of 4 or 8 channel radial catheters dependent on your preference and clinical needs.

They range from:

- Air Charged (disposable and easy to setup)
- Electronic (reusable and precise)
- Fluid Filled (reusable and low cost)

Our easy to use ARM software utilizes the same global leading software platform as our UDS software and outputs a comprehensive list of data that allows physicians to perform one of the most comprehensive ARM diagnostic tests available in the market. A variety of cabling and setup packages is available; feel free to contact your local representative for more information.

CATHETERS

Item no.	Type	Description	Recommended Uses	Quantity / Box
CAT003	Air Charged	4 Channel Radial Catheter	1	5
CAT912	Electronic	4 Channel Radial Catheter, no filling lumen	1000*	1
CAT928	Electronic	4 Channel Radial Catheter, filling lumen	1000*	1
CAT914	Electronic	8 Channel Radial Catheter, no filling lumen	1000*	1
CAT916	Electronic	8 Channel Radial Catheter, filling lumen	1000*	1
CAT006	Fluid Filled	4 Channel Radial Catheter	-**	1
CAT008	Fluid Filled	8 Channel Radial Catheter	-**	1

*According to manufacturer


**Dependant on how the catheter is handled and disinfected (Recommended useful life of 2-4 years)

ACCESSORIES


Item no.	Description	Length	Quantity / Box
TUB101	Measurement Tubing	150cm (59.1in)	25 / 50
DIS175	Graduated Cylinder (Calibration Accessory)		1
CAT752	Rectal Balloons Large (3"X1")		1



Air Charged
4 Channel Radial Catheter
CAT003



Electronic
4 Channel Radial Catheter
CAT912



Fluid Filled
4 Channel Radial Catheter
CAT006

ULTRASOUND

Needle Guides, Covers and Gel



Needle Guide Shown NDE020

Reusable Needle Guides - Endocavity and Abdominal NDE020 and NDE021

Vital component to be used as a physician's tool for guiding a needle or catheter during diagnostic ultrasound procedures in order to perform a biopsy or precise needle placement



Non-Sterile Endocavity Probe Cover

DIS370

These transrectal/transvaginal non-sterile rolled probe covers come with four rubber elastic bands. Each cover is individually wrapped and measures 2.6cm (1in) W x 3.0cm (11.8in) L.



Ultrasound Gel

ULT050 / ULT060

NuWav ultrasound gel is a standardized hypo-allergenic and is non-irritating to the most sensitive skin. Available in 250mL or 5L.

NEEDLE GUIDES

Item no.	Description	Size	Quantity / Box
NDE020	Reusable Endocavity Needle Guide used with NuWav PRB260 Probe	15cm (5.75in)	1
NDE021	Reusable Abdominal Needle Guide used with NuWav PRB240/250 Probe	6.5cm (2.5in)	1

PROBE COVERS

Item no.	Description	Length	Quantity / Box
DIS370	Non-Sterile Probe Cover, Endocavity + (4 Rubber Bands)	2.6 X 3.0cm (1X11.8in)	1

GEL

Item no.	Description	Size	Quantity / Box
ULT050	Ultrasound Gel	250mL	1
ULT060	Ultrasound Gel	5L	1

MISCELLANEOUS

Funnels, Beakers and Commodes

FUNNELS, BEAKERS

Item no.	Description	Quantity / Box
DIS173	Graduated Urine Beaker	10
CHA172	Short Funnel	1
CHA876	Long Funnel	1
STA700	Uroflow Stand	1

COMMODOES

Item no.	Description	Max Weight	Quantity / Box
CHA171	Standard Commode Chair	136.4Kg	1
CHA181	Folding Commode Chair	113.6Kg	1

PAPER

Item no.	Description	Quantity / Box
DIS155	Thermal Printer Paper Rolls	1



Folding Commode Chair

CHA181

Made from steel, this folding commode chair is height adjustable and provides a durable option as a voiding place for Uroflow or Pressure Flow Studies.

*Funnel not included



Standard Commode Chair

CHA171

This commode provides a strong, stable voiding place for Uroflow or Pressure Flow Studies. Made of PVC plastic, it is easy to clean and will not rust or chip.

*Funnel not included



Graduated Urine Beaker

DIS173

For use during Uroflow or Pressure Flow Studies, this beaker has convenient ml markers for confirmation of fluid volumes. Designed to fit perfectly on Laborie uroflow-meters.



Uroflow Stand

STA700

Made from steel, this uroflow stand is height adjustable and provides a sanitary environment for uroflow and pressure flow studies.

*Funnels not included

■ URODYNAMICS pg. 1-6

CAT206	Twin Lumen Catheter 6 FR
CAT208	Twin Lumen Catheter 8 FR
CAT307	Triple Lumen Catheter 7 FR
CAT411	Latex Rectal Balloon Catheter 10 Fr (Adult) International
CAT506	Silicone Rectal Pediatric Balloon Cath 8 FR (Pediatric) International
CAT509	Silicone Rectal Balloon Cath 10 FR (Adult) International
CAT410	Latex Rectal Balloon Cath 10 FR (Adult) USA
CAT508	Silicone Rectal Pediatric Balloon Cath 8 FR (Pediatric) USA
CAT510	Silicone Rectal Balloon Cath 10 FR (Adult) USA
CAT108	Silicone Rectal Balloon Cath 14 FR (Adult) USA
DIS102	Luer Lock Plugs
DIS103	Nova Dome
DIS203	Nova Dome & Luer Lock Plug (USA only)
DIS130	Transducer Cartridge with Luer Lock
DIS105	Catheter Adapter (USA only)
CAT103	1-Way Catheter Adapter (USA only)
DIS302	Infusion Cuff
DIS431	Stopcock - 3 - way Universal
DIS173	Graduated Urine Beaker
CAT875	Abdominal Catheter No Lumen Rectal / Vaginal 7 Fr
CAT878	Single Sensor with Lumen + Coude tip Bladder Catheter 7 Fr
CAT880	Dual Sensor Catheter with Lumen Bladder / Urethral 7 Fr
CAT885	Dual Sensor Catheter with Lumen Bladder / Urethral (Radio-opaque) 7 Fr
CAT895	Single Sensor with Lumen Bladder Catheter 7 Fr
CAT905	Single Sensor with Lumen Bladder Catheter (Radio-opaque) 7 Fr
DIS450	Air Charged Calibrator Chamber
DIS500	Stylet 10
DIS950	Velcro Leg Stabilizer Strap
TUB101	Measurement Tubing
TUB160	Clear Monitoring Line (USA only)
TUB219	Perfusion Line with 2 connections (International only)
TUB500	Pump Tubing Infusion Line
TUB501	Pump Head Tubing
TUB576	Extention Tubing without Stopcock

■ ELECTROMYOGRAPHY pg. 7

ELE335	EMG Pediatric Gel Patch Electrodes Set of 3 (USA only)
ELE425	EMG Patch Electrodes, without Lead Cables, set of 3
ELE428	EMG Patch Electrodes, with Lead Cables, set of 3
ELE200	EMG Lead Set / Shielded (Black/Red/White)
ELE365	EMG Lead Wire Non - Shielded 1 Green (USA only)
ELE370	EMG Lead Wires Non - Shielded 2 Red(USA only)
ELE350	EMG Lead Set, Non Shielded (2 Red, 1 Green)

INDEX

■ BIOFEEDBACK / STIMULATION pg. 8-10

PRB450	Vaginal Classic Probe
PRB085	Vaginal Atrophic Probe
PRB425	Vaginal Tampon Probe
PRB055	Anal Probe
ELE428	EMG Patch Electrodes, with Lead Cables, set of 3
ELE425	EMG Patch Electrodes, without Lead Cables, set of 3
ELE625	Surface Electrode / Stimulation Reusable Set of 4
DIS173	Graduated Urine Beaker
DIS203	Nova Dome with Luer Lock Plug (USA only)
DIS102	Luer Lock Plugs
DIS103	Nova Dome
PRB113	Vaginal Manometry Probe
PRB080	Vaginal Pressure Tubing
PRB108	Anal Manometry Probe
PRB078	Anal Pressure Tubing

■ ANORECTAL MANOMETRY pg. 11

CAT003	Air Charged 4 Channel Radial Catheter
CAT912	Electronic 4 Channel Radial Catheter, no filling lumen
CAT928	Electronic 4 Channel Radial Catheter, filling lumen
CAT914	Electronic 8 Channel Radial Catheter, no filling lumen
CAT916	Electronic 8 Channel Radial Catheter, filling lumen
CAT006	Fluid Filled 4 Channel Radial Catheter
CAT008	Fluid Filled 8 Channel Radial Catheter
TUB101	Measurement Tubing
DIS175	Graduated Cylinder (Calibration Accessory)
CAT752	Rectal Balloons Large (3"X1")

■ ULTRASOUND pg. 12

DIS370	Non-Sterile Probe Cover, Endocavity + (4 Rubber Bands)
ULT050	Ultrasound Gel 250mL
ULT060	Ultrasound Gel 5L

■ MISCELLANEOUS pg. 13

DIS173	Graduated Urine Beaker
CHA172	Short Funnel for CHA171 & CHA181
CHA876	Long Funnel for CHA171 & CHA181
STA700	Uroflow Stand
CHA171	Standard Commode Chair
CHA181	Folding Commode Chair
DIS155	Thermal Printer Paper Rolls

PELVIC FLOOR SYSTEMS

Aquarius TT™

The ultimate in wireless UDS technology utilizing our sleek patient tower with hi-speed EMG sampling up to 5000 samples per second, electro-mechanical height controlled transducers and a powerful processor. Aquarius TT includes the most comprehensive software package and a custom built industrial computer. Performing Video Urodynamics studies with our state of the art Smart Capture™ technology. The Aquarius TT is the right choice for the most sophisticated urologist, urogynecologist or research institution.



Triton™

Our state of the art wireless UDS system with a patient tower that can be located close to the patient while the PC or laptop can be located anywhere in the room. Hi-speed EMG sampling up to 3000 samples per second, electro-mechanical height controlled transducers and a powerful processor. A broad range of software features are standard with additional options available including an upgrade to perform video urodynamics.



Dorado KT™

The latest innovation in UDS systems by Laborie. All components are contained in an attractive, ergonomically-sound mobile cart. Optional software packages can be customized to suit your needs. Adjustable shelves and workspace trays allow for customizable configurations. The KT cart provides a high degree of in-office portability with a lockable wheelbase.



PELVIC FLOOR SYSTEMS



Delphis KT™

The Delphis KT includes a laptop computer and all components housed on the modular and mobile KT cart. Adjustable shelves and workspace trays allow for customizable configurations. The Delphis KT has a wired connection between the laptop and UDS processor, but is Bluetooth capable as an option.



Goby KT™

Urodynamics packaged in an intuitive wireless office based system. Step by Step Voice Commands walk you through each stage of a patient study. Professional One-Page Reporting at the click of the button save valuable dictation time and needless typing. The GobyKT™ is a smart choice providing flexibility and upgradeability.



Delphis™

Laborie led the industry in 2002 introducing the Delphis utilizing wireless Bluetooth technology. Compact and affordable, it is an excellent system for a clinic. The Delphis has a wireless connection between the laptop and UDS processor which allows the user to remove the laptop from the cart during patient voiding for privacy while viewing the live study.



Delphis IP™

The Delphis IP offers basic Urodynamics with software that is simple and easy to learn. An affordable, compact and portable unit mounted on an IV pole with a wireless connection to a laptop. This system is perfect for a clinic or private practice wanting to perform basic UDS studies and requiring a small footprint.

PELVIC FLOOR SYSTEMS

NuWav™ Ultrasound

A system of ultrasound probes and unique software modules for imaging bladder, kidneys, pelvic floor or the prostate. Operates independently or integrated into a Laborie urodynamics system for a multi-function, multi-modality system that can also calculate post void residual measurements.




Urostym™

The Urostym™ Behavioral Therapy System provides the most current protocols for pelvic floor retraining. Reports are quickly and easily generated providing complete patient information for analysis, treatment, and results. The addition of Uroflowmetry rounds out this device for the total picture of pelvic floor performance. Urostym can be added as an accessory to a UDS system or purchased as a stand-alone device.



Anorectal Manometry Software

As an accessory or a stand-alone system with our various hardware packages, Laborie's ARM system allows physicians to perform the latest in manometric assessments related to fecal incontinence. It provides diagnostically encompassing protocols and calculations meeting the various needs of physicians interested in Anorectal Manometry.



Urocap™

An innovative Uroflowmeter that provides accurate recordings which are easy to produce and interpret. The unit starts printing when flow is detected and stops when patient stops voiding. Automatic calculation of Uroflow Summaries. Available in a wired configuration with thermal printer or wireless to a computer.





DISPOSABLES AND ACCESSORIES CATALOGUE



HEADQUARTERS:

LABORIE MEDICAL TECHNOLOGIES Inc **International:**
6415 Northwest Drive, Unit 11, Mississauga, Ontario L4V 1X1
Phone: +1.905.612.1170
Fax: +1.905.612.9731

LABORIE MEDICAL TECHNOLOGIES Corp **USA:**
400 Avenue D, Suite 10, Williston, Vermont, 05495
Phone: 800.522.6743, +1.802.857.1300
Fax: 802.878.1122

Contact us at marketing@laborie.com
or visit us at www.laborie.com

LABORIE REPRESENTATIVE:



SALES AND SERVICE BRANCHES:

LABORIE MEDICAL TECHNOLOGIES Inc **Canada:**
2101 Boulevard Lapinière, Brossard, Québec J4W 1L7
Phone: 888.522.6743, +1.450.671.5901
Fax: +1.450.671.7182

LABORIE MEDICAL TECHNOLOGIES LTD. **Europe Distribution Center:**
Garonor Way, Royal Portbury, Bristol BS20 7XE
Phone: +44 (0)1275 372332
Fax: +44 (0)1275 375145

LABORIE MEDICAL TECHNOLOGIES N.V. **Belgium:**
BC Waasland, Industriepark-West 75, 9100 Sint-Niklaas, Belgium
Phone: +32 3 780 17 38
Fax: +32 3 780 17 39

LABORIE MEDICAL TECHNOLOGIES **France:**
24 rue Coriolis
75012 Paris, France
Phone: +33 01 43 41 34 22
Fax: +33 01 43 41 38 69

LABORIE MEDICAL TECHNOLOGIES **China:**
Room 31B, Building 3, Century Classic
No. 45, Xiao Guan Bei Li, Chao Yang District, Beijing, 100029
Phone: +86.10.8489.8891
Fax: +86.10.8489.8511

LABORIE MEDICAL TECHNOLOGIES **Brazil:**
Al.Madeira, 258 Salas 801, 805 e 806
Alphaville - 06454-010 Barueri - SP - Brazil
Phone: +55 (11) 4134 8480
Phone: +55 (11) 0 80 0771 1787
Fax: +55 (11) 4208 3837