

Pasweg 6a 3740 Bilzen Belgium

fitvibe excel pro technical data sheet



- 3D exercise images
- ITS software !!!
- Integrated exercise protocols
- Easy to use
- Start working out immediately
- Highest frequency range 20 60Hz
- Vertical vibrations
- Biggest and strongest plate
- Maximum noise prevention
- Software & Card System
- MDD medical certification (CE0197)
- Made in Germany
- 2 years warranty



Pasweg 6a 3740 Bilzen Belgium

fitvibe excel pro technical data sheet

unique selling points



Control column

- · Column separated from the vibration plate
- Pending suspension system
- · Protects the electronic parts from vibrations
- Guaranteeing a longer lifetime
- Prevent malfunctions and breakdowns

Exercise surface

- · Largest exercise surface on the market
- Allows more variations in exercises
- 5 fixation points for straps
- Strongest plate on the market
- Holds up to 250kg
- Start / Stop button





Dampening system

- 6 Absorption units
- Reducing noise
- Max prevention of vibrations to the floor
- Fine tuning of frequency between 20-60Hz



Pasweg 6a 3740 Bilzen Belgium

fitvibe excel pro technical data sheet

ITS - integrated training system software



Exercise Overview

Exercise Detail





Free Training



Muscle Group



The operating principle Clear and self-explanatory!

Free training	Start training immediately and enter your own parameters
Indications	Select the indication and enter your level, 3D-exercise dynamic pictures will guide you throughout the training protocol.
Training per muscle group	Select the zone of your body you want to strengthen up, 3D-exercise dynamic pictures will guide you trough the training protocol.
Exercise list	Choose an exercise from a list and start your training
Insert Card	Insert your personalized training card and start training

Training Indications



Exercise





Pasweg 6a 3740 Bilzen Belgium

fitvibe excel pro technical data sheet

technical specifications



Packaging:Pallet in a cardboard box. (Dimensions: 1400x 1000x 840 mm)

Vibration plate

Vertical vibration amplitude 2 & 4 mm Ergonomic handles Start / stop on platform Extra large exercise surface Frequency 20-60Hz (adjustable per 5Hz) Frequency 20-60Hz (adjustable per 1Hz) 5 strap connection points

Handle

Integrated card reader Electronics insulated from vibration Start/stop on display Start/stop by remote control

Software

Quick start Dynamic 3-D pictures Color screen Integrated training software Training & amplitude - Time card feature

Standard accessories

Foot straps (1 pair) Hand straps (1pair) Remote control Soft mat

Technical specifications

Weight 120 kg Dimensions 1050 x 990 x 1440 mm Maximum static load 250 kg (Dynamic proportions 20%) Noise level max. 71dBA Standard colors: white, grey & black Optional colors: on request Supply voltage: 230V +/- 10%, 50/60Hz 115V +/- 10%, 50/60Hz