





Donnez Uterus Manipulator

for total laparoscopic hysterectomy


Anatomical design • Easy handling and cleaning •
Short learning curve • Cost effective

Donnez* Uterus Manipulator

for total laparoscopic hysterectomy



Insulated tip



Ergonomic handle
for easy manipulation
by the assistant

Donnez* Uterus Manipulator
for total laparoscopic hysterectomy.
Ergonomic angled handle,
including cup for exposure
of anatomical structures,
completely isolated8370.601

☐ **Easy:**

- No complicated assembly or set-up prior to use
- Cleaning and sterilisation as with standard instruments

☐ **Quick:**

The learning curve for placing and using the Uterus Manipulator is short

☐ **Cost effective:**


The device is fully reusable (autoclavable)

☐ **Standard 3.5 cm cup**

☐ **Transparent cup for visualisation during initial placement**

☐ **Cup design allows exposure of anatomical structures (laparoscopic)**

☐ **Insulated cup maintains pneumoperitoneum during colpotomy**



*Jacques Donnez, M.D., Ph.D.: Department of Gynaecology, Cliniques Universitaires Saint-Luc, Louvain-en-Woluwe, Belgium

RICHARD WOLF GmbH · 75434 Knittlingen · PF 1164 · Phone +49 70 43 35-0 · Fax +49 70 43 3 53 00 · GERMANY · info@richard-wolf.com · www.richard-wolf.com

AUSTRIA:

RICHARD WOLF Austria
Ges.m.b.H.
Wilhelminenstraße 93 a
A-1160 Wien
Tel.: 01-405 51 51
Fax: 01-405 51 51 45
info@richard-wolf.at

BELGIUM:

N.V. Endoscopie
RICHARD WOLF Belgium S.A.
Industriezone Drongen, Landegemstraat 6
B-9031 Gent-Drongen
Tel.: 09.280.81.00
Fax: 09.282.92.16
endoscopy@richard-wolf.be

FRANCE:

RICHARD WOLF France S.A.R.L.
Rue Daniel Berger
Z.A.C. La Neuville
F-51100 Reims
Tél.: 03.26.87.02.89
Fax: 03.26.87.60.33
endoscopes@richardwolf.fr

UK:

RICHARD WOLF UK Ltd
Waterside Way
Wimbledon
SW17 0HB
Tel.: 020 8944 7447
Fax: 020 8944 1311
admin@richardwolf.uk.com

USA:

RICHARD WOLF
Medical Instruments Corp.
353 Corporate Woods Parkway
Vernon Hills, Illinois 60061
Tel.: 847-913 1113
Fax: 847-913 1488
sales&marketing@richardwolfusa.com