







Table of Contents

Large Bone Saw Blades		
	Compatible with DePuy® Synthes® Large Bone Oscillating Saw System	264 - 269
Accessories - Orthopaedic Pin/Drill Packs, Quick Connect Twist Drills		
	Compatible with DePuy® Synthes® Knee Instrumentation Systems	270
		
	Compatible with DePuy® Synthes® Quick Connect Couplings	271 - 273

BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K5-1270-50F	mm	70	12.5	1.27	1.27
	in	2.76	0.49	0.050	0.050
Replaces N/A					
					No Tooth Offset Flared Blade
K5-1290-47	mm	90	12.5	1.19	1.19
	in	3.54	0.49	0.047	0.047
Replaces 05.003.1045					
					No Tooth Offset
K5-1290-50	mm	90	12.5	1.27	1.27
	in	3.54	0.49	0.050	0.050
Replaces 05.003.1065					
					No Tooth Offset
K5-1290-54	mm	90	12.5	1.37	1.37
	in	3.54	0.49	0.054	0.054
Replaces 05.003.1075					
					No Tooth Offset





Compatible with DePuy® Synthes®
Large Bone Oscillating Saw System



DEPUY® SYNTHES®

BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K5-1990-39	mm	90	19	1.00	1.00
	in	3.54	0.75	0.039	0.039
Replaces 05.003.1125					
					No Tooth Offset
K5-1990-47	mm	90	19	1.19	1.19
	in	3.54	0.75	0.047	0.047
Replaces 05.003.1145					
					No Tooth Offset
K5-1990-50	mm	90	19	1.27	1.27
	in	3.54	0.75	0.050	0.050
Replaces 05.003.1165					
					No Tooth Offset
K5-1990-54	mm	90	19	1.37	1.37
	in	3.54	0.75	0.054	0.054
Replaces 05.003.1175					
					No Tooth Offset





Compatible with DePuy® Synthes®
Large Bone Oscillating Saw System



DEPUY® SYNTHES®

BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K5-2190-39F	mm	90	21	0.99	0.99
	in	3.54	0.83	0.039	0.039
Replaces 05.002.3065					
				No Tooth Offset	
				Flared Blade	
K5-2190-47F	mm	90	21	1.19	1.19
	in	3.54	0.83	0.047	0.047
Replaces 05.002.3025					
				No Tooth Offset	
				Flared Blade	
K5-2190-50F	mm	90	21	1.27	1.27
	in	3.54	0.83	0.050	0.050
Replaces 05.002.3035					
				No Tooth Offset	
				Flared Blade	
K5-2190-54F	mm	90	21	1.37	1.37
	in	3.54	0.83	0.054	0.054
Replaces 05.002.3045					
				No Tooth Offset	
				Flared Blade	





Compatible with DePuy® Synthes®
Large Bone Oscillating Saw System



DEPUY® SYNTHES®

BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K5-2590-39	mm	90	25	1.00	1.00
	in	3.54	0.98	0.039	0.039
Replaces 05.003.1325					
				No Tooth Offset	
Flared Blade					
K5-2590-47	mm	90	25	1.19	1.19
	in	3.54	0.98	0.047	0.047
Replaces 05.003.1345					
				No Tooth Offset	
K5-2590-50	mm	90	25	1.27	1.27
	in	3.54	0.98	0.050	0.050
Replaces 05.003.1365					
				No Tooth Offset	
K5-2590-54	mm	90	25	1.37	1.37
	in	3.54	0.98	0.054	0.054
Replaces 05.003.1375					
				No Tooth Offset	





Compatible with DePuy® Synthes®
Large Bone Oscillating Saw System



DEPUY® SYNTHES®

BUSA® Large Bone Saw Blades - Oscillating Replacements - DePuy® Synthes®

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM519-10	mm	50	27	0.40	0.66
	in	1.97	1.06	0.016	0.026
Replaces 519.10					
					Tooth Offset
KM519-15	mm	50	14	0.40	0.66
	in	1.97	0.55	0.016	0.026
Replaces 519.15					
					Tooth Offset
KM519-17	mm	70	27	0.64	0.81
	in	2.76	1.06	0.025	0.032
Replaces 519.17					
					Tooth Offset
KM519-21	mm	70	27	0.99	1.19
	in	2.76	1.06	0.039	0.047
Replaces 519.21					
					Tooth Offset


Compatible with DePuy® Synthes® Large Bone Oscillating Saw System




DEPUY® SYNTHES®

BUSA® Large Bone Saw Blades - Oscillating Replacements - DePuy® Synthes®

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM519-53SP	mm	75	20	0.72	0.74
	in	2.95	0.79	0.028	0.029
Replaces N/A					
					Tooth Offset 1 per Pack Non-Sterile

Special:

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-1201-84	mm	86	19.5	1.27	1.27
	in	3.39	0.77	0.050	0.050
Replaces N/A (Similar to Linvatec® 5071-201)					
					No Tooth Offset Opposing Teeth Chromium Coating

Compatible with DePuy® Synthes®
Large Bone Oscillating Saw System




DEPUY® SYNTHES®

BUSA® Orthopaedic Pin / Drill Packs - Replacements - DePuy® Synthes®

Compatible with DePuy® Synthes® Knee Instrumentation Systems

Fabricated from stainless steel.

Order Number	Product Description	Packaging
BMA-S90006S	1/8 in x 3.0 in Fixation Pin (Round End)	4 Pins/Pack (Sterile)
Replaces 2490-95-000, 2490-95		
		
BMA-S90007S	1/8 in x 3.0 in Fixation Pin (Triangular End) & Twist Drills (Round End)	4 Pins & 2 Drills/Pack (Sterile)
Replaces 86-4192		
		

Compatible with DePuy® Synthes®
Knee Instrumentation Systems

DEPUY® SYNTHES®
Replacements







Items not shown at actual size.

BUSA® Quick Connect Twist Drills

Solid - Replacements - DePuy® Synthes®

Compatible with DePuy® Synthes® Quick Connect Couplings

Fabricated from stainless steel. Features one piece construction.

Order Number		Replaces		Diameter (mm)	Total Length (mm)	Effective Length (mm)
Non-Sterile (I/pk)	Sterile (I/pk)	Non-Sterile	Sterile			
KM166-310-11	KM166-310-11S	310.11	310.11S	1.1	60	35
						
KM166-310-15	KM166-310-15S	310.15	310.15S	1.5	87	61
						
KM166-310-16	KM166-310-16S	310.16	310.16S	1.5	111	86
						
KM166-310-19	KM166-310-19S	310.19	310.19S	2.0	100	75
						
KM166-310-21	KM166-310-21S	310.21	310.21S	2.0	126	100
						
KM166-310-23	KM166-310-23S	310.23	310.23S	2.5	181	156
 <i>TiN - Titanium Nitride Coating</i>						
KM166-310-25	KM166-310-25S	310.25	310.25S	2.5	111	86
 <i>TiN - Titanium Nitride Coating</i>						
KM166-310-26	KM166-310-26S	310.26	310.26S	2.7	102	76
						
KM166-310-28	KM166-310-28S	310.28	310.28S	2.7	126	100
						








Compatible with DePuy® Synthes®
Quick Connect Couplings

DEPUY® SYNTHES®
Replacements

Items not shown at actual size.

BUSA® Quick Connect Twist Drills Solid - Replacements - DePuy® Synthes®

Compatible with DePuy® Synthes® Quick Connect Couplings

Order Number		Replaces		Diameter (mm)	Total Length (mm)	Effective Length (mm)
Non-Sterile (I/pk)	Sterile (I/pk)	Non-Sterile	Sterile			
KM166-310-29	KM166-310-29S	310.29	310.29S	3.2	196	170
						
KM166-310-31	KM166-310-31S	310.31	310.31S	3.2	147	121
						
KM166-310-35	KM166-310-35S	310.35	310.35S	3.5	112	86
						
KM166-310-37	KM166-310-37S	310.37	310.37S	3.5	196	171
						
KM166-310-44	KM166-310-44S	310.44	310.44S	4.5	147	121
						
KM166-310-48	KM166-310-48S	310.48	310.48S	4.5	195	170
						
KM166-310-530	KM166-310-530S	310.530	310.530S	2.4	100	75
						

Sterilization Parameters

Steam Sterilization Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
Pre-Vacuum (Wrapped)	270°F (132°C)	4 minutes	30 minutes
Gravity (Wrapped)	270°F (132°C)	15 minutes	30 minutes

- All cycles require a minimum of 30 minutes dry time.
- Observe AAMI and hospital guidelines for steam sterilization and the sterilizer manufacturer's written instructions for cycle parameters and load configurations.

Items not shown at actual size.






Compatible with DePuy® Synthes®
Quick Connect Couplings

DEPUY® SYNTHES®
Replacements


BUSA® Quick Connect Twist Drills Cannulated - Replacements - DePuy® Synthes®

Compatible with DePuy® Synthes® Quick Connect Couplings

Fabricated from stainless steel. Features one piece construction.

Order Number	Replaces (Non-Sterile)	Diameter (mm)	Cannulation (mm)	Total Length (mm)	Effective Length (mm)
Non-Sterile (1/pk)					
KM166-310-65	310.65	3.2	1.7	170	139
					
KM166-310-66	310.66	4.5	1.7	170	139
					
KM166-310-67	310.67	2.7	1.4	160	130
					
KM166-310-68	310.68	3.5	1.4	160	130
					
KM166-310-221	310.221	2.0	1.2	150	50
					

Note: Because of their hollow cross section, cannulated drill bits are more susceptible to breakage than solid drill bits of the same diameter. Use less axial force, avoid bending and advance the drill bit slowly to minimize the possibility of breakage.



Sterilization Parameters

Steam Sterilization Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
Pre-Vacuum (Wrapped)	270°F (132°C)	4 minutes	30 minutes
Gravity (Wrapped)	270°F (132°C)	15 minutes	30 minutes

- All cycles require a minimum of 30 minutes dry time.
- Observe AAMI and hospital guidelines for steam sterilization and the sterilizer manufacturer's written instructions for cycle parameters and load configurations.

Items not shown at actual size.

Compatible with DePuy® Synthes®
Quick Connect Couplings

DEPUY® SYNTHES®
Replacements

