

Table of Contents





















Large Bone Saw Blades		
	Compatible with Linvatec® Oscillating Saw Systems: MPower®, PowerPro® & VersiPower® Plus	199 - 210
	Compatible with Linvatec® Oscillating Saw Systems: VersiPower®, Series 4, Series 3 & Orthair	211 - 213
	Compatible with Linvatec® Reciprocating Saw Systems: MPower®, PowerPro®, VersiPower® Plus, VersiPower® Series 4, Series 3 & Orthair	214 - 215
	Compatible with Linvatec® Oscillating Saw System: Maxi-Driver™	216
		
Small Bone Saw Blades/Rasps		
	Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro I00™ & Micro E®; All BUSA® Small Bone Sagittal Saw Systems	217 - 222
	Compatible with Linvatec® Oscillating Saw Systems: MicroPower®, MicroChoice®, Micro I00™ & Micro E®; All BUSA® Small Bone Oscillating Saw Systems	223 - 226
		
	Compatible with Linvatec® Reciprocating Saw Systems: MicroPower®, MicroChoice®, Micro I00™ & Micro E®; All BUSA® Small Bone Reciprocating Saw Systems	227 - 228

Table of Contents

Small Bone Saw Blades/Rasps		
	Compatible with Linvatec® Sagittal Saw Systems: Mini-Driver™ K100 & K200 (K120 Attachment)	229 - 230
	Compatible with Linvatec® Sagittal Saw Systems: Mini-Driver™ K100 & K200 (K220 Attachment)	231
Sternum Saw Blades		
	Compatible with Linvatec® Sternum Saw Systems	232
	Compatible with Linvatec® Oscillating Saw Systems: VersiPower®, Series 4, Series 3 & Orthair; Linvatec® Micro 100™ Sagittal Saw; Linvatec® Repeat Sternotomy Saw; All BUSA® Small Bone Sagittal Saw Systems	233
	Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; All BUSA® Small Bone Sagittal Saw Systems	234
High Speed Burs		
	Compatible with High Speed Systems: Linvatec® UltraPower®, Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™	235 - 241
	Compatible with Linvatec® UltraPower® High Speed Systems	242
Standard & Oral Surgery Burs		
	Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems	243 - 257
	Compatible with Linvatec® Drill: Cebotome®	258
	Compatible with Linvatec® Drill: Osteon & Ototome	259 - 262

BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with Linvatec® Oscillating Saw Systems: MPower®, PowerPro® & VersiPower® Plus

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K7-1270-50F	mm	70	12.5	1.27	1.27
	in	2.76	0.49	0.050	0.050
Substitute for Linvatec® 5071-354					
					No Tooth Offset
Flared Blade					
K7-1290-35	mm	90	12.5	0.89	0.89
	in	3.54	0.49	0.035	0.035
Substitute for Linvatec® 5071-340					
					No Tooth Offset
K7-1290-47	mm	90	12.5	1.19	1.19
	in	3.54	0.49	0.047	0.047
Substitute for Linvatec® 5071-343					
					No Tooth Offset
K7-1290-50	mm	90	12.5	1.27	1.27
	in	3.54	0.49	0.050	0.050
Substitute for Linvatec® 5071-344					
					No Tooth Offset





Compatible with Linvatec® Oscillating Saw Systems:
MPower®, PowerPro® & VersiPower® Plus



CONMED® LINVATEC®

BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with Linvatec® Oscillating Saw Systems: MPower®, PowerPro® & VersiPower® Plus

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K7-1990-35	mm	90	19	0.89	0.89
	in	3.54	0.75	0.035	0.035
Substitute for Linvatec® 5071-300					
					No Tooth Offset
K7-1990-47	mm	90	19	1.19	1.19
	in	3.54	0.75	0.047	0.047
Substitute for Linvatec® 5071-303					
					No Tooth Offset
K7-1990-50	mm	90	19	1.27	1.27
	in	3.54	0.75	0.050	0.050
Substitute for Linvatec® 5071-304					
					No Tooth Offset
K7-1990-54	mm	90	19	1.37	1.37
	in	3.54	0.75	0.054	0.054
Substitute for Linvatec® 5071-305					
					No Tooth Offset





Compatible with Linvatec® Oscillating Saw Systems:
MPower®, PowerPro® & VersiPower® Plus



CONMED® LINVATEC®

BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with Linvatec® Oscillating Saw Systems: MPower®, PowerPro® & VersiPower® Plus

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K7-2190-39F	mm	90	21	0.99	0.99
	in	3.54	0.83	0.039	0.039
<i>Substitute for Linvatec® 5071-311</i>					
					No Tooth Offset
Flared Blade					
K7-2190-47F	mm	90	21	1.19	1.19
	in	3.54	0.83	0.047	0.047
<i>Substitute for Linvatec® 5071-313</i>					
					No Tooth Offset
Flared Blade					
K7-2190-50F	mm	90	21	1.27	1.27
	in	3.54	0.83	0.050	0.050
<i>Substitute for Linvatec® 5071-314</i>					
					No Tooth Offset
Flared Blade					
K7-2190-54F	mm	90	21	1.37	1.37
	in	3.54	0.83	0.054	0.054
<i>Substitute for Linvatec® 5071-315</i>					
					No Tooth Offset
Flared Blade					





Compatible with Linvatec® Oscillating Saw Systems:
MPower®, PowerPro® & VersiPower® Plus



CONMED® LINVATEC®

BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with Linvatec® Oscillating Saw Systems: MPower®, PowerPro® & VersiPower® Plus

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K7-2590-35	mm	90	25	0.89	0.89
	in	3.54	.98	0.035	0.035
Substitute for Linvatec® 5071-320					
					No Tooth Offset
K7-2590-39	mm	90	25	0.99	0.99
	in	3.54	.98	0.039	0.039
Substitute for Linvatec® 5071-321					
					No Tooth Offset
K7-2590-47	mm	90	25	1.19	1.19
	in	3.54	.98	0.047	0.047
Substitute for Linvatec® 5071-323					
					No Tooth Offset
K7-2590-50	mm	90	25	1.27	1.27
	in	3.54	.98	0.050	0.050
Substitute for Linvatec® 5071-324					
					No Tooth Offset



Compatible with Linvatec® Oscillating Saw Systems:
MPower®, PowerPro® & VersiPower® Plus



CONMED® LINVATEC®

BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with Linvatec® Oscillating Saw Systems: MPower®, PowerPro® & VersiPower® Plus

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K7-2590-54	mm	90	25	1.37	1.37
	in	3.54	.98	0.054	0.054
Substitute for Linvatec® 5071-325					
					No Tooth Offset
K7-2590-58	mm	90	25	1.47	1.47
	in	3.54	.98	0.058	0.058
Substitute for Linvatec® 5071-326					
					No Tooth Offset





Compatible with Linvatec® Oscillating Saw Systems:
MPower®, PowerPro® & VersiPower® Plus



CONMED® LINVATEC®

BUSA® Large Bone Saw Blades - Oscillating - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: MPower®, PowerPro® & VersiPower® Plus

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM71-201	mm	88	19.5	1.27	1.27
	in	3.44	0.77	0.050	0.050
Replaces 5071-201					
			No Tooth Offset		
			Opposing Teeth Chromium Coating		
KM71-233	mm	25	9.5	0.40	0.60
	in	.98	0.37	0.016	0.024
Replaces 5071-133, 5071-233					
			Tooth Offset		
			ACL Blade		
KM71-523	mm	97	19.5	1.27	1.27
	in	3.79	0.77	0.050	0.050
Replaces 5071-123, 5071-523					
			No Tooth Offset		
			Opposing Teeth		
KM71-524	mm	97	25.4	1.27	1.27
	in	3.79	1.00	0.050	0.050
Replaces 5071-124, 5071-524					
			No Tooth Offset		
			Opposing Teeth		





Compatible with Linvatec® Oscillating Saw Systems:
MPower®, PowerPro® & VersiPower® Plus



CONMED® LINVATEC®
Replacements

BUSA® Large Bone Saw Blades - Oscillating - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: MPower®, PowerPro® & VersiPower® Plus

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM71-552	mm	65	13.5	0.61	0.90
	in	2.56	0.53	0.024	0.035
Replaces 5071-152, 5071-452, 5071-552					
					Tooth Offset Flared Blade
KM71-553	mm	65	19.5	0.61	0.90
	in	2.56	0.77	0.024	0.035
Replaces 5071-153, 5071-453, 5071-553					
					Tooth Offset Opposing Teeth
KM71-554	mm	70	19.5	1.20	1.40
	in	2.76	0.77	0.047	0.055
Replaces 5071-154, 5071-454, 5071-554					
					Tooth Offset Opposing Teeth
KM71-556	mm	65	33	0.80	0.90
	in	2.56	1.30	0.031	0.035
Replaces 5071-1456, 5071-456, 5071-556					
					Tooth Offset Opposing Teeth





Compatible with Linvatec® Oscillating Saw Systems:
MPower®, PowerPro® & VersiPower® Plus



CONMED® LINVATEC®
Replacements

BUSA® Large Bone Saw Blades - Oscillating - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: MPower®, PowerPro® & VersiPower® Plus

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM71-561	mm	92	25.4	1.27	1.27
	in	3.59	1.00	0.050	0.050
Replaces 5071-161, 5071-561					
					No Tooth Offset
					Opposing Teeth
KM71-563	mm	92	25.4	1.37	1.37
	in	3.59	1.00	0.054	0.054
Replaces 5071-163, 5071-563					
					No Tooth Offset
					Opposing Teeth
KM71-567	mm	92	19.5	1.20	1.20
	in	3.59	0.77	0.047	0.047
Replaces 5071-167, 5071-567					
					No Tooth Offset
					Opposing Teeth Chromium Coating
KM71-572	mm	70	19.5	0.80	0.90
	in	2.76	0.77	0.031	0.035
Replaces 5071-172, 5071-472, 5071-572					
					Tooth Offset
					Opposing Teeth





Compatible with Linvatec® Oscillating Saw Systems:
MPower®, PowerPro® & VersiPower® Plus



CONMED® LINVATEC®
Replacements

BUSA® Large Bone Saw Blades - Oscillating - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: MPower®, PowerPro® & VersiPower® Plus

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM71-574	mm	85	19.5	0.80	0.90
	in	3.35	0.77	0.031	0.035
Replaces 5071-174, 5071-574					
					Tooth Offset Opposing Teeth
KM71-577	mm	70	19.5	0.80	0.90
	in	2.76	0.77	0.031	0.035
Replaces 5071-177, 5071-577					
					Tooth Offset Opposing Teeth Chromium Coating
KM71-578	mm	80	19.5	0.80	0.90
	in	3.15	0.77	0.031	0.035
Replaces 5071-178, 5071-578					
					Tooth Offset Opposing Teeth Chromium Coating
KM71-580	mm	85	10	1.27	1.27
	in	3.35	0.39	0.050	0.050
Replaces 5071-180, 5071-580					
					No Tooth Offset Diamond Tooth Design





Compatible with Linvatec® Oscillating Saw Systems:
MPower®, PowerPro® & VersiPower® Plus



CONMED® LINVATEC®
Replacements

BUSA® Large Bone Saw Blades - Oscillating - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: MPower®, PowerPro® & VersiPower® Plus

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM71-581	mm	88	19.5	1.27	1.27
	in	3.44	0.77	0.050	0.050
Replaces 5071-181, 5071-581					
					No Tooth Offset
					Opposing Teeth
KM71-585	mm	82	19.5	0.90	0.90
	in	3.20	0.77	0.035	0.035
Replaces 5071-185, 5071-585					
					No Tooth Offset
					Opposing Teeth
KM71-586	mm	65	32	0.80	0.90
	in	2.56	1.26	0.031	0.035
Replaces 5071-186, 5071-586					
					Tooth Offset
					Opposing Teeth
KM71-595	mm	82	19.5	0.90	0.90
	in	3.20	0.77	0.035	0.035
Replaces 5071-195, 5071-595					
					No Tooth Offset
					Opposing Teeth Chromium Coating


Compatible with Linvatec® Oscillating Saw Systems:
MPower®, PowerPro® & VersiPower® Plus



CONMED® LINVATEC®
Replacements

BUSA® Large Bone Saw Blades - Oscillating - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: MPower®, PowerPro® & VersiPower® Plus

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM71-596	mm	92	25.4	0.90	0.90
	in	3.59	1.00	0.035	0.035
Replaces 5071-196, 5071-596					
 No Tooth Offset Opposing Teeth					





Compatible with Linvatec® Oscillating Saw Systems:
MPower®, PowerPro® & VersiPower® Plus



CONMED® LINVATEC®
Replacements

BUSA® Large Bone Saw Blades - Oscillating - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: VersiPower®, Series 4, Series 3 & Orthair

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-251	mm	60	32.5	0.61	0.90
	in	2.36	1.28	0.024	0.035
Replaces 5052-251, 5052-451, 5052-551					
					Tooth Offset Flared Blade Opposing Teeth
KM-256	mm	60	33.5	0.80	0.90
	in	2.36	1.32	0.031	0.035
Replaces 5052-256, 5052-456, 5052-556					
					Tooth Offset Opposing Teeth
KM-272	mm	70	19	0.80	0.90
	in	2.76	0.76	0.031	0.035
Replaces 5052-272, 5052-472, 5052-572					
					Tooth Offset Opposing Teeth
KM-274	mm	85	19	0.80	0.90
	in	3.35	0.76	0.031	0.035
Replaces 5052-274, 5052-474, 5052-574					
					Tooth Offset Opposing Teeth

Compatible with Linvatec® Oscillating Saw Systems:
VersiPower®, Series 4, Series 3 & Orthair



CONMED® LINVATEC®
Replacements





BUSA® Large Bone Saw Blades - Oscillating - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: VersiPower®, Series 4, Series 3 & Orthair

Compatible with Linvatec® Oscillating Saw Systems:
VersiPower®, Series 4, Series 3 & Orthair





CONMED® LINVATEC®
Replacements

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-450	mm	65	49.1	0.61	0.90
	in	2.56	1.93	0.024	0.035
Replaces 5052-250, 5052-450, 5052-550					
					Tooth Offset Flared Blade
KM-451	mm	60	32.5	0.61	0.90
	in	2.36	1.28	0.024	0.035
Replaces 5052-251, 5052-451, 5052-551					
					Tooth Offset Flared Blade
KM-452	mm	60	12.5	0.61	0.90
	in	2.36	0.49	0.024	0.035
Replaces 5052-252, 5052-452, 5052-552					
					Tooth Offset Flared Blade
KM-453	mm	60	19.5	0.61	0.92
	in	2.36	0.77	0.024	0.036
Replaces 5052-253, 5052-453, 5052-553					
					Tooth Offset Opposing Teeth

BUSA® Large Bone Saw Blades - Oscillating - Replacements - Linvatec®


Compatible with Linvatec® Oscillating Saw Systems: VersiPower®, Series 4, Series 3 & Orthair

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-454	mm	70	19.5	1.27	1.50
	in	2.76	0.77	0.050	0.059
Replaces 5052-254, 5052-454, 5052-554					
					Tooth Offset
KM-1201-44	mm	88	19.5	1.27	1.27
	in	3.44	0.77	0.050	0.050
Replaces N/A (Similar to Linvatec® 5071-201)					
					No Tooth Offset Opposing Teeth Chromium Coating

Compatible with Linvatec® Oscillating Saw Systems:
VersiPower®, Series 4, Series 3 & Orthair

BUSA® Large Bone Saw Blades - Oscillating - Replacements - Synvasive®

Compatible with Linvatec® Oscillating Saw Systems: VersiPower®, Series 4, Series 3 & Orthair




Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM11-0030	mm	90	25.4	1.27	1.27
	in	3.54	1.00	0.050	0.050
Replaces 11-0030					
					No Tooth Offset Opposing Teeth



CONMED® LINVATEC®
Replacements



BUSA® Large Bone Saw Blades-Reciprocating -Replacements- Linvatec®

Compatible with Linvatec® Reciprocating Saw Systems: MPower®, PowerPro®, VersiPower® Plus, VersiPower® Series 4, Series 3 & Orthair

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-458	mm	89	12.5	0.61	1.00
	in	3.50	0.49	0.024	0.039
Replaces 5052-058, 5052-258					
					Tooth Offset
KM-459	mm	54.1	9.4	0.61	1.00
	in	2.13	0.37	0.024	0.039
Replaces 5052-059, 5052-259					
					Tooth Offset
KM-460	mm	38	9.6	0.40	0.60
	in	1.50	0.38	0.016	0.024
Replaces 5052-060, 5052-260					
					Tooth Offset
KM-461	mm	89	12.5	1.50	1.60
	in	3.50	0.49	0.059	0.063
Replaces 5052-061, 5052-261					
					Tooth Offset

BUSA® Large Bone Saw Blades-Reciprocating -Replacements- Linvatec®

Compatible with Linvatec® Reciprocating Saw Systems: MPower®, PowerPro®, VersiPower® Plus, VersiPower® Series 4, Series 3 & Orthair





Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-473	mm	73	12.5	0.80	1.00
	in	2.87	0.49	0.031	0.039
Replaces 5052-053, 5052-173					
					Tooth Offset
					Double Sided
KM-479	mm	77.5	12.7	0.90	0.90
	in	3.05	0.50	0.035	0.035
Replaces 5052-079, 5052-179					
					Tooth Offset
					Double Sided
					Ribbed
					Chromium Coating

Compatible with Linvatec® Reciprocating Saw Systems: MPower®, PowerPro®, VersiPower® Plus, VersiPower® Series 4, Series 3 & Orthair

CONMED® LINVATEC®
Replacements

BUSA® Large Bone Saw Blades - Oscillating - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw System: Maxi-Driver™





Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-22L	mm	60	29	0.97	1.62
	in	2.36	1.16	0.038	0.064
Replaces L122					
					Tooth Offset
KM-23L	mm	60	14.7	0.61	1.17
	in	2.36	0.59	0.024	0.046
Replaces L123					
					Tooth Offset
KM-24L	mm	55	31	0.61	1.14
	in	2.17	1.22	0.024	0.045
Replaces L124					
					Tooth Offset
KM-27L	mm	55	31	0.97	1.49
	in	2.17	1.22	0.038	0.059
Replaces L127					
					Tooth Offset

Compatible with Linvatec® Oscillating
Saw System: Maxi-Driver™

CONMED® LINVATEC®
Replacements

BUSA® Small Bone Sagittal Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro I 00™ & Micro E®;
All BUSA® Small Bone Sagittal Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3100C	mm	41.5	14	0.40	0.61
	in	1.63	0.55	0.016	0.024
Replaces 5023-133, 5023-233 Replaces 5053-033, 5053-233					
					Coarse Straight Edge
KM-3101	mm	25.7	9.5	0.40	0.60
	in	1.01	0.37	0.016	0.024
Replaces 5023-138, 5023-238 Replaces 5053-038, 5053-238					
					Fine Straight Edge
KM-3101C	mm	25.7	9.5	0.40	0.61
	in	1.01	0.37	0.016	0.024
Replaces 5023-143, 5023-243 Replaces 5053-043, 5053-243					
					Coarse Straight Edge
KM-3102	mm	25.7	4.5	0.40	0.60
	in	1.01	0.18	0.016	0.024
Replaces 5023-139, 5023-239 Replaces 5053-039, 5053-239					
					Fine Straight Edge





Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro I 00™ & Micro E®, All BUSA® Small Bone Sagittal Saw Systems



CONMED® LINVATEC®
Replacements





BUSA® Small Bone Sagittal Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro I 00™ & Micro E®; All BUSA® Small Bone Sagittal Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3104	mm	16.1	9.5	0.40	0.63
	in	0.63	0.37	0.016	0.025
Replaces N/A					
 Fine Straight Edge					
KM-3105	mm	16.1	4.5	0.40	0.63
	in	0.63	0.18	0.016	0.025
Replaces N/A					
 Fine Straight Edge					
KM-3106	mm	25.7	5.5	0.40	0.61
	in	1.01	0.22	0.016	0.024
Replaces 5023-132, 5023-232 Replaces 5053-032, 5053-232					
 Coarse Straight Edge					
KM-3107	mm	34.3	6	0.40	0.60
	in	1.35	0.24	0.016	0.024
Replaces 5023-141, 5023-241 Replaces 5053-041, 5053-241					
 Fine Straight Edge					

BUSA® Small Bone Sagittal Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®;
All BUSA® Small Bone Sagittal Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3183	mm	25.7	9.5	0.40	0.60
	in	1.01	0.37	0.016	0.024
Replaces 5023-183					
					Fine Straight Edge ACL Blade
KM-3231	mm	41.5	19.5	0.40	0.61
	in	1.63	0.77	0.016	0.024
Replaces 5023-131, 5023-231 Replaces 5053-031, 5053-231					
					Coarse Straight Edge
KM-3237	mm	25.7	14	0.40	0.61
	in	1.01	0.55	0.016	0.024
Replaces 5023-137, 5023-237 Replaces 5053-037, 5053-237					
					Coarse Straight Edge
KM-3240	mm	41.5	9.5	0.40	0.60
	in	1.63	0.37	0.016	0.024
Replaces 5023-140, 5023-240 Replaces 5053-040, 5053-240					
					Fine Straight Edge Tapered



Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®, All BUSA® Small Bone Sagittal Saw Systems



CONMED® LINVATEC®
Replacements






BUSA® Small Bone Sagittal Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro I 00™ & Micro E®;
All BUSA® Small Bone Sagittal Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3179SN	mm	42.5	54	0.40	0.69
	in	1.67	2.13	0.016	0.027
Replaces N/A					
 Coarse Opposing Teeth Non-Sterile					
KM-3231SN	mm	41.9	19	0.40	0.71
	in	1.65	0.75	0.016	0.028
Replaces N/A					
 Non-Sterile					

BUSA® Small Bone Sagittal Saw Blades Graft Harvesting - Replacements - Linvatec®

Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®;
All BUSA® Small Bone Sagittal Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3405	mm	20	16.5	0.40	5.0
	in	0.79	0.65	0.016	0.197
Replaces 5023-405					
  Graft Harvesting Blade					
KM-3406	mm	20	16.5	0.40	6.0
	in	0.79	0.65	0.016	0.236
Replaces 5023-406					
  Graft Harvesting Blade					
KM-3407	mm	20	16.5	0.40	7.0
	in	0.79	0.65	0.016	0.276
Replaces 5023-407					
  Graft Harvesting Blade					
KM-3408	mm	20	16.5	0.40	8.0
	in	0.79	0.65	0.016	0.315
Replaces 5023-408					
  Graft Harvesting Blade					

Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®, All BUSA® Small Bone Sagittal Saw Systems







CONMED® LINVATEC®
Replacements







BUSA® Small Bone Sagittal Saw Blades Graft Harvesting - Replacements - Linvatec®

Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro I 00™ & Micro E®;
All BUSA® Small Bone Sagittal Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3409	mm	20	16.5	0.40	9.0
	in	0.79	0.65	0.016	0.354
Replaces 5023-409					
  Graft Harvesting Blade					
KM-3410	mm	20	16.5	0.40	10.0
	in	0.79	0.65	0.016	0.394
Replaces 5023-410					
  Graft Harvesting Blade					

BUSA® Small Bone Oscillating Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®;
All BUSA® Small Bone Oscillating Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3000	mm	18.5	5.5	0.40	0.60
	in	0.73	0.22	0.016	0.024
Replaces 5023-162, 5023-262 Replaces 5053-062, 5053-262					
 Fine Straight Edge					
KM-3001	mm	28.5	5.5	0.40	0.60
	in	1.12	0.22	0.016	0.024
Replaces 5023-164, 5023-264 Replaces 5053-064, 5053-264					
 Fine Straight Edge					
KM-3002	mm	34.7	10	0.40	0.60
	in	1.37	0.39	0.016	0.024
Replaces 5023-169, 5023-269 Replaces 5053-073, 5053-273					
 Fine Straight Edge					
KM-3003	mm	22	5.5	0.40	0.60
	in	0.87	0.22	0.016	0.024
Replaces N/A					
 Fine Straight Edge					





Compatible with Linvatec® Oscillating Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; All BUSA® Small Bone Oscillating Saw Systems



CONMED® LINVATEC®
Replacements



BUSA® Small Bone Oscillating Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; All BUSA® Small Bone Oscillating Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3004	mm	27.9	9.5	0.40	0.60
	in	1.10	0.37	0.016	0.024
Replaces N/A					
 Fine Straight Edge					
KM-3005	mm	18.5	10	0.40	0.60
	in	0.73	0.39	0.016	0.024
Replaces 5023-163, 5023-263 Replaces 5053-063, 5053-263					
 Fine Straight Edge					
KM-3006	mm	12	5.5	0.40	0.60
	in	0.47	0.22	0.016	0.024
Replaces 5023-160, 5023-260 Replaces 5053-060, 5053-260					
 Fine Straight Edge					
KM-3008	mm	28.5	18.5	0.40	0.60
	in	1.12	0.73	0.016	0.024
Replaces 5023-166, 5023-266 Replaces 5053-066, 5053-266					
 Fine Radial Edge					

BUSA® Small Bone Oscillating Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®;
All BUSA® Small Bone Oscillating Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3261	mm	12	10.0	0.40	0.60
	in	0.47	0.39	0.016	0.024
Replaces 5023-161, 5023-261 Replaces 5053-061, 5053-261					
 Fine Straight Edge					
KM-3265	mm	28.5	10.00	0.40	0.60
	in	1.12	0.39	0.016	0.024
Replaces 5023-165, 5023-265 Replaces 5053-065, 5053-265					
 Fine Radial Edge					




Compatible with Linvatec® Oscillating Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; All BUSA® Small Bone Oscillating Saw Systems



CONMED® LINVATEC®
Replacements



BUSA® Small Bone Oscillating Saw Blades Crescentic - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: MicroPower®, MicroChoice®, Micro I 00™ & Micro E®;
All BUSA® Small Bone Oscillating Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3170	mm	12	9	0.40	0.60
	in	0.47	0.36	0.016	0.024
Replaces 5023-170, 5023-270 Replaces 5053-170, 5053-270					
 Fine Crescentic Blade					
KM-3171	mm	25.5	13.8	0.40	0.60
	in	1.00	0.55	0.016	0.024
Replaces 5023-171, 5023-271 Replaces 5053-171, 5053-271					
 Fine Crescentic Blade					
KM-3172	mm	32	18	0.40	0.60
	in	1.26	0.71	0.016	0.024
Replaces 5023-172, 5023-272 Replaces 5053-172, 5053-272					
 Fine Crescentic Blade					
KM-3176	mm	32	13.8	0.40	0.60
	in	1.26	0.55	0.016	0.024
Replaces 5023-176, 5023-276 Replaces 5053-176, 5053-276					
 Fine Crescentic Blade					

BUSA® Small Bone Reciprocating Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Reciprocating Saw Systems: MicroPower®, MicroChoice®, Micro I 00™ & Micro E®;
All BUSA® Small Bone Reciprocating Saw Systems





Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-060	mm	25.5	7.6	0.40	0.64
	in	1.00	0.30	0.016	0.025
Replaces 5053-020, 5053-220					
					Flat Shank
KM-070	mm	38	7.6	0.40	0.64
	in	1.50	0.30	0.016	0.025
Replaces 5053-021, 5053-221					
					Flat Shank
KM-54	mm	23.2	7.4	0.40	0.60
	in	0.913	0.29	0.016	0.024
Replaces 5053-054, 5053-254					
					Round Shank

Compatible with Linvatec® Reciprocating Saw Systems: MicroPower®, MicroChoice®, Micro I 00™ & Micro E®, All BUSA® Small Bone Reciprocating Saw Systems

CONMED® LINVATEC®
Replacements






BUSA® Small Bone Reciprocating Rasps - Replacements - Linvatec®

Compatible with Linvatec® Reciprocating Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; All BUSA® Small Bone Reciprocating Saw Systems

Order Number	Unit	Length	Max. Width
KM-80	mm	11	5.5
	in	0.43	0.22
Replaces 5053-080			
 Small Tear Round Shank			
KM-81	mm	13	6.9
	in	0.51	0.27
Replaces 5053-081			
 Large Tear Round Shank			
KM-82	mm	11	5.5
	in	0.43	0.22
Replaces 5053-082			
 Small Tear Cross-Cut Round Shank			
KM-83	mm	13	6.9
	in	0.51	0.27
Replaces 5053-083			
 Large Tear Cross-Cut Round Shank			

BUSA® Small Bone Sagittal Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Sagittal Saw Systems: Mini-Driver™ K100 & K200 (K120 Attachment)

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-21K	mm	15.1	6.4	0.71	0.89
	in	0.59	0.25	0.028	0.035
Replaces K121					
 Fine					
KM-22K	mm	31	6.4	0.71	0.89
	in	1.22	0.25	0.028	0.035
Replaces K122					
 Fine					
KM-23K	mm	15.1	10.4	0.71	0.89
	in	0.59	0.41	0.028	0.035
Replaces K123					
 Fine					
KM-24K	mm	31	10.4	0.71	0.89
	in	1.22	0.41	0.028	0.035
Replaces K124					
 Fine					
KM-31K	mm	15.1	6.4	0.38	0.56
	in	0.59	0.25	0.015	0.022
Replaces K131					
 Fine					



Compatible with Linvatec® Sagittal Saw Systems:
Mini-Driver™ K100 & K200 (K120 Attachment)



CONMED® LINVATEC®
Replacements

BUSA® Small Bone Sagittal Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Sagittal Saw Systems: Mini-Driver™ K100 & K200 (K120 Attachment)

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-33K	mm	15.1	10.4	0.38	0.56
	in	0.59	0.41	0.015	0.022
Replaces K133					
 Fine					
KM-34K	mm	31	10.4	0.38	0.56
	in	1.22	0.41	0.015	0.022
Replaces K134					
 Fine					






Compatible with Linvatec® Sagittal Saw Systems:
Mini-Driver™ K100 & K200 (K120 Attachment)



CONMED® LINVATEC®
Replacements

BUSA® Small Bone Sagittal Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Sagittal Saw Systems: Mini-Driver™ K100 & K200 (K220 Attachment)



Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-41K	mm	14.2	6.4	0.71	0.97
	in	0.56	0.25	0.028	0.038
Replaces K141					
 Fine					
KM-43K	mm	14.2	10.2	0.71	0.97
	in	0.56	0.40	0.028	0.038
Replaces K143					
 Fine					
KM-44K	mm	30	10.2	0.71	0.97
	in	1.18	0.40	0.028	0.038
Replaces K144					
 Fine					
KM-48K	mm	27.4	9.5	0.38	0.65
	in	1.08	0.37	0.015	0.026
Replaces K148					
 Fine					
KM-51K	mm	14.2	6.4	0.38	0.65
	in	0.56	0.25	0.015	0.026
Replaces K151					
 Fine					

Compatible with Linvatec® Sagittal Saw Systems:
Mini-Driver™ K100 & K200 (K220 Attachment)

 **CONMED® LINVATEC®**
Replacements

BUSA® Sternum Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Sternum Saw Systems



Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-31	mm	54.0	9.9	0.64	.97
	in	2.13	0.39	0.025	0.038
Replaces 5059-31, 5059-031, 5059-531, 5059-031-01					
					Sterile
KM-31N	mm	54.0	9.9	0.64	.97
	in	2.13	0.39	0.025	0.038
Replaces 5059-31, 5059-031, 5059-531, 5059-031-01					
					Non-Sterile
KM-32	mm	32.8	9.4	0.64	.99
	in	1.29	0.37	0.025	0.039
Replaces 5059-32, 5059-032, 5059-532, 5059-032-01					
					Sterile
KM-32N	mm	32.8	9.4	0.64	.99
	in	1.29	0.37	0.025	0.039
Replaces 5059-32, 5059-032, 5059-532, 5059-032-01					
					Non-Sterile

Compatible with Linvatec®
Sternum Saw Systems

CONMED® LINVATEC®
Replacements

BUSA® Sternum Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: VersiPower®, Series 4, Series 3 & Orthair; Linvatec® Micro 100™ Sagittal Saw; Linvatec® Repeat Sternotomy Saw; All BUSA® Small Bone Sagittal Saw Systems



Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-279	mm	39.5	55.0	0.40	0.79
Replaces 5059-279	in	1.55	2.15	0.016	0.031
 Sterile Repeat Sternotomy Blade Opposing Teeth					
KM-279PZ	mm	26.0	39.0	0.40	0.71
Replaces N/A	in	1.02	1.53	0.016	0.028
 Sterile Repeat Sternotomy Blade					

Compatible with Linvatec® Oscillating Saw Systems: VersiPower®, Series 4, Series 3 & Orthair; Linvatec® Micro 100™ Sagittal Saw; Linvatec® Repeat Sternotomy Saw; All BUSA® Small Bone Sagittal Saw Systems


CONMED® LINVATEC®
 Replacements

BUSA® Sternum Saw Blades - Replacements - Linvatec®


Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; All BUSA® Small Bone Sagittal Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3179SN	mm	42.5	54	0.40	0.69
	in	1.67	2.13	0.016	0.027
Replaces N/A					
 Coarse Opposing Teeth Non-Sterile					
KM-3231SN	mm	41.9	19	0.40	0.71
	in	1.65	0.75	0.016	0.028
Replaces N/A					
 Non-Sterile					

BUSA® High Speed Burs - Replacements

Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-AM2U	Acorn	mm	7.5	9.4
		in	.30	.37
Replaces Linvatec® 7021-382 Replaces Medtronic® Midas Rex® AM-2 Replaces Stryker® 5130-30-75, 5230-30-75, 5430-30-75				
				
KM-M8U	Match Head	mm	2.9	3.8
		in	0.11	0.15
Replaces Linvatec® 7021-237 Replaces Medtronic® Midas Rex® M-8 Replaces Stryker® 5120-107-30, 5220-107-30, 5420-107-30				
				
KM-AM8U	Match Head	mm	2.9	3.8
		in	0.11	0.15
Replaces Linvatec® 7021-337 Replaces Medtronic® Midas Rex® AM-8 Replaces Stryker® 5130-107-30, 5230-107-30, 5430-107-30				
				
KM-AM8DU	Match Head Diamond	mm	2.9	3.8
		in	0.11	0.15
Replaces Linvatec® 7021-379 Replaces Medtronic® Midas Rex® AM-8D Replaces Stryker® 5130-107-130, 5230-107-130, 5430-107-32				
				

Compatible with High Speed Systems: Linvatec® UltraPower®, Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™






CONMED® LINVATEC®
Replacements



BUSA® High Speed Burs - Replacements Linvatec®, Medtronic® Midas Rex® & Stryker®


Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-M51U	Metal Cutter (Diamond Wheel)	mm	25.4	0.8
		in	1.00	0.03
Replaces Linvatec® 7021-292 Replaces Medtronic® Midas Rex® M-51 Replaces Stryker® 5120-130-250				
				
KM-M52U	Metal Cutter (Carbide)	mm	3.1	17.5
		in	0.12	0.69
Replaces Linvatec® 7021-294 Replaces Medtronic® Midas Rex® M-52 Replaces Stryker® 5120-80-30, 5420-80-30				
				
KM-U295	Oval	mm	4.0	8.0
		in	0.16	0.31
Replaces Linvatec® 7021-295 Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5120-15-140				
				
KM5110-15-40	Oval	mm	4.0	9.5
		in	0.16	0.37
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5110-15-40				
				

BUSA® High Speed Burs - Replacements

Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM5130-15-40	Oval	mm	4.0	9.5
		in	0.16	0.37
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5130-15-40				
				
KM-M340DU	Round Diamond	mm	4.0	4.0
		in	0.16	0.16
Replaces Linvatec® 7021-264 Replaces Medtronic® Midas Rex® M-340D Replaces Stryker® 5120-12-40				
				
KM-U265	Round Diamond	mm	5.0	5.0
		in	0.20	0.20
Replaces Linvatec® 7021-265 Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5120-12-50				
				
KM-M320U	Round Fluted	mm	2.0	2.0
		in	0.08	0.08
Replaces Linvatec® 7021-224 Replaces Medtronic® Midas Rex® M-320 Replaces Stryker® 5120-10-20, 5220-10-20				
				

Compatible with High Speed Systems: Linvatec® UltraPower®, Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™





CONMED® LINVATEC®
Replacements





BUSA® High Speed Burs - Replacements Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-M330U	Round Fluted	mm	3.0	3.0
		in	0.12	0.12
Replaces Linvatec® 7021-226 Replaces Medtronic® Midas Rex® M-330 Replaces Stryker® 5120-10-30, 5220-10-30				
				
KM5120-10-230	Round Fluted Extended Length	mm	3.0	3.0
		in	0.12	0.12
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5120-10-230, 5220-10-230				
				
KM-M340U	Round Fluted	mm	4.0	4.0
		in	0.16	0.16
Replaces Linvatec® 7021-227 Replaces Medtronic® Midas Rex® M-340 Replaces Stryker® 5120-10-40, 5220-10-40				
				
KM5120-10-240	Round Fluted Extended Length	mm	4.0	4.0
		in	0.16	0.16
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5120-10-240, 5220-10-240				
				

BUSA® High Speed Burs - Replacements Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-M35U	Round Fluted	mm	5.0	5.0
		in	0.20	0.20
Replaces Linvatec® 7021-228 Replaces Medtronic® Midas Rex® M-35 Replaces Stryker® 5120-10-50, 5220-10-50, 5420-10-50				
				
KM5120-10-250	Round Fluted Extended Length	mm	5.0	5.0
		in	0.20	0.20
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5120-10-250, 5220-10-250				
				
KM-M33U	Round Fluted	mm	6.0	6.0
		in	0.24	0.24
Replaces Linvatec® 7021-229 Replaces Medtronic® Midas Rex® M-33 Replaces Stryker® 5120-10-60, 5220-10-60				
				
KM-M32U	Round Fluted	mm	7.6	7.6
		in	0.30	0.30
Replaces Linvatec® 7021-231 Replaces Medtronic® Midas Rex® M-32 Replaces Stryker® 5120-10-75, 5220-10-75				
				





Compatible with High Speed Systems: Linvatec® UltraPower®, Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™



CONMED® LINVATEC®
Replacements

BUSA® High Speed Burs - Replacements Linvatec®, Medtronic® Midas Rex® & Stryker®




Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-M31U	Round Fluted	mm	9.0	9.0
		in	0.35	0.35
Replaces Linvatec® 7021-232 Replaces Medtronic® Midas Rex® M-31 Replaces Stryker® 5120-10-90, 5220-10-90				
				
KM-AM330U	Round Fluted	mm	3.0	3.0
		in	0.12	0.12
Replaces Linvatec® 7021-326 Replaces Medtronic® Midas Rex® AM-330 Replaces Stryker® 5130-10-30, 5230-10-30				
				
KM-AM340U	Round Fluted	mm	4.0	4.0
		in	0.16	0.16
Replaces Linvatec® 7021-327 Replaces Medtronic® Midas Rex® AM-340 Replaces Stryker® 5130-10-40, 5230-10-40				
				
KM5130-10-240	Round Fluted Extended Length	mm	4.0	4.0
		in	0.16	0.16
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5130-10-240, 5230-10-240				
				

BUSA® High Speed Burs - Replacements

Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-AM35U	Round Fluted	mm	5.0	5.0
		in	0.20	0.20
Replaces Linvatec® 7021-328 Replaces Medtronic® Midas Rex® AM-35 Replaces Stryker® 5130-10-50, 5230-10-50, 5430-10-50				
				
KM5130-10-250	Round Fluted Extended Length	mm	5.0	5.0
		in	0.20	0.20
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5130-10-250, 5230-10-250				
				
KM-AM33U	Round Fluted	mm	6.0	6.0
		in	0.24	0.24
Replaces Linvatec® 7021-329 Replaces Medtronic® Midas Rex® AM-33 Replaces Stryker® 5130-10-60, 5230-10-60				
				
KM-AM32U	Round Fluted	mm	7.6	7.6
		in	0.30	0.30
Replaces Linvatec® 7021-331 Replaces Medtronic® Midas Rex® AM-32 Replaces Stryker® 5130-10-75, 5230-10-75				
				

Compatible with High Speed Systems: Linvatec® UltraPower®, Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™



CONMED® LINVATEC®
Replacements

BUSA® High Speed Burs - Replacements - Linvatec®

Compatible with Linvatec® UltraPower® High Speed Systems

Order Number	Description	Unit	Head Diameter	Head Length
KM-U815	Craniotome Bur	mm	2.2	16.0
		in	0.09	0.63
Replaces Linvatec® 7021-815 Replaces Medtronic® Midas Rex® N/A Replaces Stryker® N/A				
				

Compatible with Linvatec®
UltraPower® High Speed Systems







CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery

2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems





Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®, MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-C1803	Lindemann - Carbide	mm	51.0	1.9	13.0	4
		in	2.01	0.07	0.51	
Replaces N/A						
						
03-C0901	Oval - Carbide	mm	48.0	4.0	8.0	8
		in	1.89	0.16	0.31	
Replaces 091-136, 5091-236						
						
03-C0905	Oval - Carbide	mm	48.0	5.5	10.2	10
		in	1.89	0.22	0.40	
Replaces 5091-138, 5091-238						
						
09-C0901	Oval - Carbide	mm	70.0	4.0	8.0	8
		in	2.76	0.16	0.31	
Replaces 5092-136, 5092-236						
						
09-C0905	Oval - Carbide	mm	70.0	5.5	10.2	10
		in	2.76	0.22	0.40	
Replaces 5092-138, 5092-238						
						
12-C0901	Oval - Carbide	mm	95.0	4.0	8.0	8
		in	3.74	0.16	0.31	
Replaces 5093-136, 5093-236						
						

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems


Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
12-C0905	Oval - Carbide	mm	95.0	5.5	10.2	10
		in	3.74	0.22	0.40	
Replaces 5093-138, 5093-238						
						
04-C0101	Round - Carbide	mm	51.0	1	1	6
		in	2.01	0.04	0.04	
Replaces 5091-120, 5091-220						
						
04-C0103	Round - Carbide	mm	51.0	1.5	1.5	6
		in	2.01	0.06	0.06	
Replaces 5091-122, 5091-222						
						
04-C0104	Round - Carbide	mm	51.0	1.8	1.8	6
		in	2.01	0.07	0.07	
Replaces 5091-121, 5091-221						
						
04-C0105	Round - Carbide	mm	51.0	2.0	2.0	6
		in	2.01	0.08	0.08	
Replaces 5091-124, 5091-224						
						
04-C0106	Round - Carbide	mm	51.0	2.3	2.3	8
		in	2.01	0.09	0.09	
Replaces 5091-123, 5091-125, 5091-223, 5091-225						
						

BUSA® Micro Burs - Standard & Oral Surgery

2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems



Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®, MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-C0107	Round - Carbide	mm	51.0	3.0	3.0	8
		in	2.01	0.12	0.12	
Replaces 5091-126, 5091-226						
						
04-C0108	Round - Carbide	mm	51.0	4.0	4.0	8
		in	2.01	0.16	0.16	
Replaces 5091-128, 5091-228						
						
04-C0109	Round - Carbide	mm	51.0	5.0	5.0	10
		in	2.01	0.20	0.20	
Replaces 5091-130, 5091-230						
						
04-C0110	Round - Carbide	mm	51.0	6.0	6.0	8
		in	2.01	0.24	0.24	
Replaces M60R						
						
04-C0170	Round - Carbide	mm	51.0	6.5	6.5	10
		in	2.01	0.26	0.26	
Replaces 5091-132, 5091-232						
						
09-C0109	Round - Carbide	mm	70.0	5.0	5.0	10
		in	2.76	0.20	0.20	
Replaces 5092-130, 5092-230						
						

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems






Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
09-C0110	Round - Carbide	mm	70.0	6.0	6.0	8
		in	2.76	0.24	0.24	
Replaces L60R						
						
09-C0170	Round - Carbide	mm	70.0	6.5	6.5	10
		in	2.76	0.26	0.26	
Replaces 5092-132, 5092-232						
						
10-C0101	Round - Carbide	mm	71.5	1.0	1.0	6
		in	2.81	0.04	0.04	
Replaces 5092-120, 5092-220						
						
10-C0103	Round - Carbide	mm	72.0	1.5	1.5	6
		in	2.83	0.06	0.06	
Replaces 5092-122, 5092-222						
						
10-C0105	Round - Carbide	mm	72.5	2.0	2.0	6
		in	2.85	0.08	0.08	
Replaces 5092-124, 5092-224						
						
10-C0107	Round - Carbide	mm	73.0	3.0	3.0	8
		in	2.87	0.12	0.12	
Replaces 5092-126, 5092-226						
						

BUSA® Micro Burs - Standard & Oral Surgery

2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems





Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®, MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
10-C0108	Round - Carbide	mm	73.0	4.0	4.0	8
		in	2.87	0.16	0.16	
Replaces 5092-128, 5092-228						
						
12-C0105	Round - Carbide	mm	95.0	2.0	2.0	6
		in	3.74	0.08	0.08	
Replaces 5093-124, 5093-224						
						
12-C0108	Round - Carbide	mm	95.0	4.0	4.0	8
		in	3.74	0.16	0.16	
Replaces 5093-128, 5093-228						
						
12-C0109	Round - Carbide	mm	95.0	5.0	5.0	10
		in	3.74	0.20	0.20	
Replaces 5093-130, 5093-230						
						
12-C0170	Round - Carbide	mm	95.0	6.5	6.5	10
		in	3.74	0.26	0.26	
Replaces 5093-132, 5093-232						
						
04-C0601	Tapered - Carbide	mm	51.0	0.9	5.2	6
		in	2.01	0.035	0.20	
Replaces 5091-115, 5091-215						
						

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems



Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-C0604	Tapered - Carbide	mm	51.0	1.0	6.0	6
		in	2.01	0.04	0.24	
Replaces N/A						
						
04-C0606	Tapered - Carbide	mm	51.0	1.2	4.2	6
		in	2.01	0.05	0.17	
Replaces N/A						
						
04-C0623	Tapered - Carbide	mm	51.0	1.2	6.0	6
		in	2.01	0.05	0.24	
Replaces N/A						
						
04-C0609	Tapered - Carbide	mm	51.0	1.5	10.0	6
		in	2.01	0.06	0.39	
Replaces 5091-102, 5091-202						
						
04-C0607	Tapered - Carbide	mm	51.0	1.6	4.3	6
		in	2.01	0.06	0.17	
Replaces 5091-118, 5091-218						
						
04-C0617	Tapered - Carbide	mm	51.0	2.0	7.5	6
		in	2.01	0.08	0.30	
Replaces 5091-106, 5091-206						
						

BUSA® Micro Burs - Standard & Oral Surgery

2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems



Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®, MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-C0615	Tapered - Carbide	mm	51.0	2.1	4.8	6
		in	2.01	0.08	0.19	
Replaces N/A						
						
04-C0612	Tapered - Carbide	mm	51.0	2.1	4.9	6
		in	2.01	0.08	0.19	
Replaces 5091-113, 5091-213						
						
06-C0619	Tapered - Carbide	mm	57.0	2.0	19.0	6
		in	2.24	0.08	0.75	
Replaces 5091-114, 5091-214						
						
09-C0613	Tapered - Carbide	mm	70.0	1.73	12.0	6
		in	2.76	0.07	0.47	
Replaces 5092-103, 5092-203						
						
10-C0617	Tapered - Carbide	mm	73.0	2.0	7.5	6
		in	2.87	0.08	0.30	
Replaces 5092-106, 5092-206						
						
11-C0619	Tapered - Carbide	mm	80.0	2.0	19.0	6
		in	3.15	0.08	0.75	
Replaces 5092-114, 5092-214						
						

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems






Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
03-C0312	Straight - Cross-Cut - Carbide	mm	48.0	4.0	10.0	6
		in	1.89	0.16	0.39	
Replaces 5091-154, 5091-254						
						
04-C0302	Straight - Cross-Cut - Carbide	mm	51.0	1.0	5.0	6
		in	2.01	0.04	0.20	
Replaces 5091-101, 5091-201						
						
04-C0404	Straight - Cross-Cut - Carbide	mm	51.0	1.2	3.7	6
		in	2.01	0.05	0.15	
Replaces 5091-119, 5091-219						
						
04-C0306	Straight - Cross-Cut - Carbide	mm	51.0	1.5	7.5	6
		in	2.01	0.06	0.30	
Replaces 5091-104, 5091-204						
						
04-C0307	Straight - Carbide	mm	51.0	1.6	4.3	6
		in	2.01	0.06	0.17	
Replaces 5091-117, 5091-217						
						
04-C0310	Straight - Cross-Cut - Carbide	mm	51.0	2.04	10.0	8
		in	2.01	0.08	0.39	
Replaces 5091-152, 5091-252						
						

BUSA® Micro Burs - Standard & Oral Surgery

2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems






Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®, MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
10-C0306	Straight - Carbide	mm	73.0	1.5	7.5	6
		in	2.87	0.06	0.30	
Replaces 5092-104, 5092-204						
						
10-C0310	Straight - Cross-Cut - Carbide	mm	73.0	2.04	10.0	8
		in	2.87	0.08	0.39	
Replaces 5092-152, 5092-252						
						
03-D0911*	Oval - Diamond - Dermabrader	mm	48.0	4.0	8.0	Medium
		in	1.89	0.16	0.31	
Replaces 5091-170, 5091-270						
						
*Check Manufacturer's recommended bur guards and speed for handpiece being used.						
03-D0135	Round - Diamond	mm	48.0	4.0	4.0	Medium
		in	1.89	0.16	0.16	
Replaces 5091-164, 5091-264						
						
03-D0138	Round - Diamond	mm	48.0	6.0	6.0	Coarse
		in	1.89	0.24	0.24	
Replaces 5091-168, 5091-268						
						
04-D0139	Round - Diamond	mm	51.0	2.0	2.0	Medium
		in	2.01	0.08	0.08	
Replaces 5091-160, 5091-260						
						

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems


Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-D0133	Round - Diamond	mm	51.0	3.0	3.0	Medium
		in	2.01	0.12	0.12	
Replaces 5091-162, 5091-262						
						
09-D0128	Round - Diamond	mm	70.0	1.0	1.0	Fine
		in	2.76	0.04	0.04	
Replaces N/A						
						
09-D0140	Round - Diamond	mm	70.0	2.3	2.3	Medium
		in	2.76	0.09	0.09	
Replaces N/A						
						
09-D0132	Round - Diamond	mm	70.0	2.9	2.9	Medium
		in	2.76	0.11	0.11	
Replaces N/A						
						
09-D0135	Round - Diamond	mm	70.0	4.0	4.0	Medium
		in	2.76	0.16	0.16	
Replaces 5092-164, 5092-264						
						
09-D0141	Round - Diamond	mm	70.0	5.0	5.0	Coarse
		in	2.76	0.20	0.20	
Replaces 5092-166, 5092-266						
						

BUSA® Micro Burs - Standard & Oral Surgery

2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems



Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®, MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
09-D0138	Round - Diamond	mm	70.0	6.0	6.0	Coarse
		in	2.76	0.24	0.24	
Replaces 5092-168, 5092-268						
						
10-D0139	Round - Diamond	mm	73.0	2.0	2.0	Medium
		in	2.87	0.08	0.08	
Replaces 5092-160, 5092-260						
						
10-D0133	Round - Diamond	mm	73.0	3.0	3.0	Medium
		in	2.87	0.12	0.12	
Replaces 5092-162, 5092-262						
						
12-D0133	Round - Diamond	mm	95.0	3.0	3.0	Medium
		in	3.74	0.12	0.12	
Replaces 5093-162, 5093-262						
						
12-D0135	Round - Diamond	mm	95.0	4.0	4.0	Medium
		in	3.74	0.16	0.16	
Replaces 5093-164, 5093-264						
						
12-D0141	Round - Diamond	mm	95.0	5.0	5.0	Coarse
		in	3.74	0.20	0.20	
Replaces 5093-166, 5093-266						
						

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems





Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
08-S1802	Lindemann - Stainless Steel	mm	67.0	2.1	25.0	4
		in	2.64	0.083	.98	
Replaces 5091-112, 5091-212						
						
03-S1101	Flame - Stainless Steel	mm	48.0	6.0	13.0	14
		in	1.89	0.24	0.51	
Replaces N/A						
						
09-S0101	Round - Stainless Steel	mm	70.0	1.0	1.0	6
		in	2.76	0.04	0.04	
Replaces N/A						
						
09-S0150	Round - Stainless Steel	mm	70.0	1.4	1.4	8
		in	2.76	0.058	0.058	
Replaces N/A						
						
09-S0151	Round - Stainless Steel	mm	70.0	1.8	1.8	8
		in	2.76	0.07	0.07	
Replaces N/A						
						
09-S0115	Round - Stainless Steel	mm	70.0	2.4	2.4	8
		in	2.76	0.091	0.091	
Replaces N/A						
						

BUSA® Micro Burs - Standard & Oral Surgery

2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems





Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®, MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
09-S0153	Round - Stainless Steel	mm	70.0	3.1	3.1	10
		in	2.76	0.12	0.12	
Replaces N/A						
						
09-S0120	Round - Stainless Steel	mm	70.0	4.0	4.0	12
		in	2.76	0.16	0.16	
Replaces N/A						
						
09-S0155	Round - Stainless Steel	mm	70.0	5.0	5.0	12
		in	2.76	0.20	0.20	
Replaces N/A						
						
09-S0156	Round - Stainless Steel	mm	70.0	6.0	6.0	14
		in	2.76	0.24	0.24	
Replaces N/A						
						
09-S0125	Round - Stainless Steel	mm	70.0	7.0	7.0	16
		in	2.76	0.28	0.28	
Replaces N/A						
						
12-S0153	Round - Stainless Steel	mm	95.0	3.1	3.1	10
		in	3.74	0.12	0.12	
Replaces N/A						
						

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems


Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
12-S0123	Round - Stainless Steel	mm	95.0	4.5	4.5	12
		in	3.74	0.179	0.18	
Replaces N/A						
						
12-S0155	Round - Stainless Steel	mm	95.0	5.0	5.0	12
		in	3.74	0.20	0.20	
Replaces N/A						
						
04-S1501	Shannon Style - Tapered - Stainless Steel	mm	51.0	2.3	12.0	4
		in	2.01	0.09	0.47	
Replaces N/A						
						
04-S1502	Shannon Style - Straight - Stainless Steel	mm	51.0	2.4	13.2	3
		in	2.01	0.091	0.52	
Replaces N/A						
						
06-S1601	Swanson Reamer - Stainless Steel	mm	57.0	2.0	10.0	5
		in	2.24	0.08	0.39	
Replaces 5091-142, 5091-242						
						
07-S1602	Swanson Reamer - Stainless Steel	mm	62.0	3.0	13.5	5
		in	2.44	0.12	0.53	
Replaces 5091-143, 5091-243						
						

BUSA® Micro Burs - Standard & Oral Surgery

2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems



Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®, MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
07-S1605	Swanson Reamer - Stainless Steel	mm	62.0	4.0	16.5	6
		in	2.44	0.16	0.65	
Replaces 5091-144, 5091-244						
						
01-S1703	Wire Pass Drill - Stainless Steel	mm	37.0	1.5	19.0	2
		in	1.46	0.059	0.75	
Replaces 5090-148, 5090-248						
						
01-S1704	Wire Pass Drill - Stainless Steel	mm	37.0	2.0	19.0	2
		in	1.46	0.079	0.75	
Replaces 5090-149, 5090-249						
						
04-S1701	Wire Pass Drill - Stainless Steel	mm	51.0	1.0	19.0	2
		in	2.01	0.039	0.75	
Replaces 5091-147, 5091-247						
						
04-S1703	Wire Pass Drill - Stainless Steel	mm	51.0	1.5	19.0	2
		in	2.01	0.059	0.75	
Replaces 5091-148, 5091-248						
						
04-S1704	Wire Pass Drill - Stainless Steel	mm	51.0	2.0	19.0	2
		in	2.01	0.079	0.75	
Replaces 5091-149, 5091-249						
						

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard Cebotome® Shank - Replacements - Linvatec®

Compatible with Linvatec® Drill: Cebotome®



Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
13-S01189	Round - Stainless Steel	mm	30.1	5.0	5.0	6
		in	1.19	0.20	0.20	
Replaces 5052-189						
						
13-S01184	Round - Stainless Steel	mm	30.3	7.0	7.0	10
		in	1.19	0.28	0.28	
Replaces 5052-184						
						
13-S01185	Round - Stainless Steel	mm	31.5	8.0	8.0	10
		in	1.24	0.31	0.31	
Replaces 5052-185						
						
13-S01186	Round - Stainless Steel	mm	32.6	9.0	9.0	10
		in	1.28	0.35	0.35	
Replaces 5052-186						
						

Compatible with Linvatec®
Drill: Cebotome®

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard Osteon/Ototome Shank - Replacements - Linvatec®

Compatible with Linvatec® Drill: Osteon & Ototome







Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
70-C0108	Round - Carbide	mm	30.0	4.0	4.0	8
		in	1.18	0.16	0.16	
Replaces 5056-19, 5056-119						
						
70-C0109	Round - Carbide	mm	30.0	5.0	5.0	10
		in	1.18	0.20	0.20	
Replaces 5056-18, 5056-118						
						
70-C0170	Round - Carbide	mm	30.0	6.5	6.5	10
		in	1.18	0.26	0.26	
Replaces 5056-17, 5056-117						
						
71-C0107	Round - Carbide	mm	38.0	3.0	3.0	8
		in	1.50	0.12	0.12	
Replaces 5056-28, 5056-128						
						
71-C0108	Round - Carbide	mm	38.0	4.0	4.0	8
		in	1.50	0.16	0.16	
Replaces 5056-27, 5056-127						
						
71-C0109	Round - Carbide	mm	38.0	5.0	5.0	10
		in	1.50	0.20	0.20	
Replaces 5056-26, 5056-126						
						

Compatible with Linvatec®
Drill: Osteon & Ototome

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard Osteon/Ototome Shank - Replacements - Linvatec®

Compatible with Linvatec® Drill: Osteon & Ototome







Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
71-C0170	Round - Carbide	mm	38.0	6.5	6.5	10
		in	1.50	0.26	0.26	
Replaces 5056-25, 5056-125						
						
72-C0107	Round - Carbide	mm	48.0	3.0	3.0	8
		in	1.89	0.12	0.12	
Replaces 5056-46, 5056-146						
						
72-C0108	Round - Carbide	mm	48.0	4.0	4.0	8
		in	1.89	0.16	0.16	
Replaces 5056-45, 5056-145						
						
72-C0109	Round - Carbide	mm	48.0	5.0	5.0	10
		in	1.89	0.20	0.20	
Replaces 5056-44, 5056-144						
						
75-S0101	Round - Stainless Steel	mm	48.0	1.0	1.0	6
		in	1.89	0.04	0.04	
Replaces 5056-64, 5056-164						
						
75-S0148	Round - Stainless Steel	mm	48.0	0.56	0.5	6
		in	1.89	0.022	0.02	
Replaces 5056-66, 5056-166						
						

Compatible with Linvatec®
Drill: Osteon & Ototome

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard Osteon/Ototome Shank - Replacements - Linvatec®

Compatible with Linvatec® Drill: Osteon & Ototome





Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
75-S0149	Round - Stainless Steel	mm	48.0	0.74	0.74	6
		in	1.89	0.03	0.03	
Replaces 5056-65, 5056-165						
						
75-S0150	Round - Stainless Steel	mm	48.0	1.4	1.4	8
		in	1.89	0.058	0.058	
Replaces 5056-40, 5056-140						
						
75-S0151	Round - Stainless Steel	mm	48.0	1.8	1.8	8
		in	1.89	0.07	0.07	
Replaces 5056-42, 5056-142						
						
71-D0135	Round - Diamond	mm	38.0	4.0	4.0	Medium
		in	1.50	0.16	0.16	
Replaces 5056-48, 5056-148						
						
71-D0141	Round - Diamond	mm	38.0	5.0	5.0	Coarse
		in	1.50	0.20	0.20	
Replaces 5056-50, 5056-150						
						
72-D0133	Round - Diamond	mm	48.0	3.0	3.0	Medium
		in	1.89	0.12	0.12	
Replaces 5056-59, 5056-159						
						

Compatible with Linvatec®
Drill: Osteon & Ototome

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard Osteon/Ototome Shank - Replacements - Linvatec®

Compatible with Linvatec® Drill: Osteon & Ototome

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
72-D0135	Round - Diamond	mm	48.0	4.0	4.0	Medium
		in	1.89	0.16	0.16	
Replaces 5056-52, 5056-152						
						
72-D0141	Round - Diamond	mm	48.0	5.0	5.0	Coarse
		in	1.89	0.20	0.20	
Replaces N/A						
						
75-D0128	Round - Diamond	mm	48.0	1.0	1.0	Fine
		in	1.89	0.04	0.04	
Replaces 5056-58, 5056-158						
						
75-D0139	Round - Diamond	mm	48.0	2.0	2.0	Medium
		in	1.89	0.08	0.08	
Replaces 5056-62, 5056-162						
						
75-D0142	Round - Diamond	mm	48.0	0.54	0.54	X-Fine
		in	1.89	0.021	0.021	
Replaces 5056-61, 5056-161						
						
75-D0143	Round - Diamond	mm	48.0	0.72	0.72	Fine
		in	1.89	0.03	0.03	
Replaces 5056-60, 5056-160						
						

Compatible with Linvatec®
Drill: Osteon & Ototome

CONMED® LINVATEC®
Replacements