



IV Pole for infusion pumps and infusion fluids

Link to product:

<http://meden.com.pl/offer/covid-19/228-iv-pole-for-infusion-pumps-and-infusion-fluids.html>

Condition: New

Informacje:

IV Pole is a functional stand for hanging containers with infusion fluid. A thoughtful design and a stable 5-arm base (made of aluminum or high-quality plastic) facilitate everyday use and movement.

Product description:

IV Pole for infusion pumps and infusion fluids

Stand is distinguished by high durability, reliability, mobility and ease of use. And is great for various types of medical facilities (hospitals, clinics) and doctor's offices and for home.

The fluids hanger can be locked at the desired height using the knob located on the column. On the column can be easily mounted infusion pumps (the types which are equipped with a universal holder for vertical stands).

Features:

- Solid and durable construction made of powder-coated stainless steel,
- Hanger with 2 holders for infusion fluids,
- Easy height adjustment with a knob,
- Possibility of mounting several infusion pumps,
- 5-arm base made of aluminum or plastic,
- Rubber coated wheels for easy and quiet movement.





Technical data:

Height [cm]:	136-240
Width [cm]:	52
Weight [kg]:	6
Max. load of the holders (load of 1 holder / load of 2 holders) [kg]:	3/8
Max. load of the stand column [kg]:	10
Max. safe stand load (SWL) [kg]:	18