



Adjustable hygiene chair Mohican

Link to product: <http://meden.com.pl/offer/home/150-adjustable-hygiene-chair-mohican.html>

Available online only

Condition: New

Informacje:

Mobile hygiene chair with **adjustable height** and positioning is protected against water, dirt or hair while bathing, shaving or grooming the patient. Adjusted with a **manual remote controller** for optimal positioning of the patient.

[View product presentation](#)



Product description:

Mohican allows for an easy access to the **shower** and helps in **moving the patient around** at home or in care facilities of all kinds. Suitable **for the elderly** or **people with all kinds of disabilities** who might need help in taking care of their **personal hygiene**.

Product features:

- grips for easy movement,
- waterproof, easy to clean backrest,

- PU foam seat with an opening to help with patient hygiene,
 - two waterproofed, dustproofed, battery-powered actuators,
 - manual remote controller,
 - 4 blockable wheels,
 - optional headrest for maximum comfort of the patient.
-

Optional accessories:

- adjustable chest belt - easy to use,
 - velcro straps for legs - easy to remove and clean,
 - extra battery - high-capacity pack for longer use.
-

Technical data:

| | |
|--|-----------|
| Size (LxW) [cm]: | 110 x 57 |
| Distance between armrests [cm]: | 47 |
| Seat tilt [°]: | +7 to -30 |
| Max. load [kg]: | 150 |
| Weight [kg]: | 44 |
| Height adjustment [cm]: | 53-93 |
| Power [V]: | 24 |