



## Adjustable hygiene chair Mohican

Link to product: <http://meden.com.pl/offer/home/150-adjustable-hygiene-chair-mohican.html>

Available online only

Condition: New

### Informacje:

**Mobile hygiene chair** with **adjustable height** and positioning is protected against water, dirt or hair while bathing, shaving or grooming the patient. Adjusted with a **manual remote controller** for optimal positioning of the patient.

[View product presentation](#)



### Product description:

**Mohican** allows for an easy access to the **shower** and helps in **moving the patient around** at home or in care facilities of all kinds. Suitable **for the elderly** or **people with all kinds of disabilities** who might need help in taking care of their **personal hygiene**.

### Product features:

- grips for easy movement,
- waterproof, easy to clean backrest,

- PU foam seat with an opening to help with patient hygiene,
  - two waterproofed, dustproofed, battery-powered actuators,
  - manual remote controller,
  - 4 blockable wheels,
  - optional headrest for maximum comfort of the patient.
- 

**Optional accessories:**

- adjustable chest belt - easy to use,
  - velcro straps for legs - easy to remove and clean,
  - extra battery - high-capacity pack for longer use.
- 

**Technical data:**

<b>Size (LxW) [cm]:</b>	110 x 57
<b>Distance between armrests [cm]:</b>	47
<b>Seat tilt [°]:</b>	+7 to -30
<b>Max. load [kg]:</b>	150
<b>Weight [kg]:</b>	44
<b>Height adjustment [cm]:</b>	53-93
<b>Power [V]:</b>	24