



Hydromassage membrane unit Aquai

Link to product:

<http://meden.com.pl/offer/hydrotherapy-medical-equipment/1-hydromassage-membrane-unit-aquai.html>

Available online only

Condition: New

Informacje:

Dry hydromassage membrane unit. The **AQUAI** has wide range of uses: it can be applied in various **wellness** and **fitness** areas. It perfectly suits for **rehabilitation centers**. It is also a great proposal for **relaxation rooms** in companies.

[View product presentation](#)



Product description:

Hydromassage membrane unit Aquai

Advantages of the treatment: improvement in **blood circulation** and **metabolism** increase in **oxygen transport** in muscle enhancement of the production of **endorphins** responsible for the **pain tolerance** improvement in **movement range** of extremities **tension relief stress reduction** and **relaxation well-being**

Treatment possibilities:

A patient in a **rehabilitation center** as well as a guest in a **SPA center** can experience beneficial power of hydrojet massage **without getting wet**. This effect is possible thanks to the layer of **plastic coat** through which the massage is done. The user is placed on the **elastic membrane** which bottom side is hit by the water **running from the nozzles**. Regardless of the number of treatments, the **Aquai** is able to maintain constant water temperature. For the additional comfort, the **Aquai** can memorize a specific **massage program** the user wishes to use more frequently.

Advantages of the treatment:

- improvement in blood circulation and metabolism
- increase in oxygen transport in muscle
- enhancement of the production of endorphins responsible for the pain tolerance
- improvement in movement range of extremities
- tension relief
- stress reduction and relaxation
- well-being

Product features:

- 6 programmable and independent zones
- all jets are able to work simultaneously
- 6 standard massage programs
- 1 user program
- membrane surface temperature setting: temperature (adjustable from 30 to 40°C)
- can be programmed for determined time
- adjustable treatment time (1 - 30 min)

- delayed start at programmed number or hours
 - adjustable frequency and sequence of zones
 - adjustable water temperature
 - integrated heating system
 - no need to connect to a water tap
 - easily accessible touch panel
 - LED lighting at the base
-

Technical data:

Capacity:	220 l
Power consumption:	~3,5 kW
Supply voltage:	230 V / 50 Hz or 400V / 50-60 Hz
Dimensions:	240 cm x 115 cm x 65 cm
Weight:	190 kg empty/410 kg filled with water

Watch product presentation: