



Jet shower ORKAN type NBC1

Link to product:

<http://meden.com.pl/offer/hydrotherapy-medical-equipment/71-jet-shower-orkan-type-nbc1.html>

Available online only

Condition: New

Informacje:

Jet shower with two water jet guns.

Application:

- treatment of muscle and joint rheumatism while there is no inflammation
- treatment of sciatica
- treatment of stomach and bowel disfunctions



Product description:

Orkan whip showers allow to use thermal and mechanical water properties during therapeutic treatment. Using stimuli provided by the stream of water and rapid water temperature changes, the therapeutic procedure can affect patient in different ways, for example: stimulate or sedate, increase or decrease the blood pressure.

Product features:

- Free-standing console
- 2 manometers and 2 thermometers (independent for each jet gun),
- Water jet dispersion adjustment (independent for both jet guns),
- The ability to work with two jet guns at the same time,
- The jet guns rotatable within the treatment area,
- The jet guns equipped with triggers initiating water massage (with a lock allowing for continuous massage) with forearm support,
- 2 thermostatic taps for independent temperature adjustment for the two jet guns,
- 2 water massage intensity (pressure) adjustment and shut-off valves independent for the two jet guns.

Technical data:

Warm and cold water pressure [bar]:	4-6
Warm water temperature [°C]:	max. 60
Power supply [V/Hz]:	230/50
Installed power [kW]:	0,9
Dimensions NBC1 (H x W x L) [cm]:	115 x 67 x 55