

KTM, KTM BO

Fine motoric tables for manual hand exercises



Box dimensions: **80 x 30 x 60 cm** Weight: **+1,5 kg**

KTM:

- The exercises on the tables improve the ability to daily activities
- Therapeutic rehabilitation in case of orthopedic, rheumatological and neurological disease of hand
- For use in the hospitals (in rehabilitation, surgical or neurological wards) as well as in outpatient rehabilitation
- Improve the manipulation skills, motor coordination of hand and eye-hand coordination
- For dorsiflexion, palmar flexion and rotation in wrist joint, conversion and inversion of the forearm
- Foldable base

Accessories:

- Circle with handle - for conversion and inversion of the forearm
- Handle - for dorsiflexion and palmar flexion in wrist joint
- Flap - for dorsiflexion and palmar flexion in wrist joint
- Vertical spiral - for rotation in wrist joint
- Horizontal roller - for dorsiflexion and palmar flexion in wrist joint
- Arm stabilizer (stabilizer axis adjustment in the range of 0-90°)

Technical data:

Table top dimensions (W x L) [cm]:	52 x 86
Height adjustment [cm]:	54-86
Mass [kg]:	10
Plummet mass (5 items) [kg]:	0,25/item

KTM BO:

- The exercises on the tables improve the ability to daily activities
- Therapeutic rehabilitation in case of orthopedic, rheumatological and neurological disease of hand
- For use in the hospitals (in rehabilitation, surgical or neurological wards) as well as in outpatient rehabilitation
- Improve the manipulation skills, motor coordination of hand and eye-hand coordination
- For dorsiflexion, palmar flexion and rotation in wrist joint, conversion and inversion of the forearm
- Foldable base

Accessories:

- Horizontal spiral - for rotation in wrist joint
- Vertical roller with the spring - for dorsiflexion and palmar flexion in wrist joint
- Vertical ball with the spring - for abduction and adduction in wrist joint
- Pin - for improvement of manipulatin skills
- Knob - for improvement of manipulatin skills
- Arm stabilizer (stabilizer axis adjustment in the range of 0-90°)

Technical data:

Table top dimensions (W x L) [cm]:	72 x 52
Height adjustment [cm]:	54-86
Mass [kg]:	8,4

