



## Solmed UVA-B Ultraviolet irradiation lamp

Link to product:

<http://meden.com.pl/offer/light-therapy-equipment-phototherapy/168-solmed-uva-b-ultraviolet-irradiation-lamp.html>

Available online only

Condition: New

### Informacje:

**Solmed UVA-B** is an ultraviolet irradiation lamp intended for treatment with the use of wide spectrum of ultraviolet rays.

### Product description:

## Solmed UVA-B Ultraviolet irradiation lamp

The lamp is applicable to the treatment of certain affections, improves metabolism, improves bone remineralisation, stimulates endocrine glands, improves oxygen metabolism, lowers hypertension and has beneficial influence on the haematopoiesis. Intended users are physiotherapy, rehabilitation units and sanatoriums.

### Functions:

- Irradiation
- Adjustable treatment time
- Adjustment of the radiators position
- Adjustment of the radiators height

### Main indications to the use of UV light:

- Skin disorders (acne vulgaris etc.)
- Furuncles
- Inflammation of soft tissue
- Alopecia areata
- Thropic ulcerations, slow-healing wounds
- Psoriasis
- Respiratory system conditions (nose and throat conditions, chronic bronchitis, asthma), laryngological diseases
- Rheumatoid diseases (soft tissue rheumatic conditions, rheumatoid arthritis, arthrosis)
- During convalescence
- Hypofunction of endocrine glands
- Hindered bone tissue adhesion
- Prevention of rickets

### Technical data:

Duration of treatment [min.]:	1-15 (option 1-60)
Power supply [V/Hz]:	220-240/50-60
Max. single radiator power [W]:	20
Max. power consumption[A]:	0,82
Fuses:	2 x T2AH/250 V
Dimensions (L x W x H) [cm]:	62 x 48 x 187
Weight [kg]:	32