



Phototherapy device Solmed UNO

Link to product:

<http://meden.com.pl/offer/light-therapy-equipment-phototherapy/49-phototherapy-device-solmed-uno.html>

Available online only

Condition: New

Informacje:

Irradiation lamp with 1 radiator is intended for treatment by heating the skin tissue using the heat produced by the infrared light. The lamp can be used to cure some of the chronic and subacute inflammatory states and pain syndromes. It can also be used as preparation treatment before a massage, manual therapy, iontophoresis or kinesitherapy.

Product description:

Product features:

- infrared light emission
- 1 radiator in the head of the unit
- adjustable head
- height of the stand adjusted by means of a gas spring
- light source equipped with a cooling fan and a filter mounting system
- variable radiation power
- adjustable duration of the treatment
- control panel with buttons: low-profile wheeled base for easy moving, which fits under massage tables
 - START /PAUSE
 - STOP
 - SETUP button (adjusting the power and duration of treatment)
 - PROGRAM button (programming the sequence of treatment by setting up different power levels for each time interval)

Standard components:

- unit with head (1)
- bulb (1)
- red filter (1)
- blue filter (1)
- big protective goggles for the operator (1)
- small protective goggles for the patient (1)
- assembly kit
- instruction manual

Technical data:

Power supply:	230 V / 50 Hz
Maximum bulb power:	375 W
Power consumption:	max. 395 W
Dimensions:	740 x 480 x 1760 mm
Weight:	about 26 kg
Duration of treatment:	30 min
Power adjustment:	10 - 100 %

