



Fine motoric tables without resistance for manual hand exercises KTM BO

Link to product:


<http://meden.com.pl/offer/limb-dysfunction-therapy/149-fine-motoric-tables-without-resistance-f-or-manual-hand-exercises-ktm-bo.html>

Available online only

Condition: New

Informacje:

Fine motoric tables without resistance for manual hand exercises.



Product description:

Fine motoric tables without resistance for manual hand exercises.

Application:

- the exercises on the tables improve the ability to daily activities,
- therapeutic rehabilitation in case of orthopedic, rheumatological and neurological disease of hand,
- for use in the hospitals (in rehabilitation, surgical or neurological wards) as well as in outpatient rehabilitation,
- improve the manipulation skills, motor coordination of hand and eye-hand coordination,
- for dorsiflexion, palmar flexion and rotation in wrist joint, conversion and inversion of the forearm.

Characteristics:

- simply tables - both for adults and children (height adjustment 54 - 86 cm),
- the design of the table is made of wood and steel,
- top painted with ecological varnish,
- 4 non-slip rubber feet with height adjustment,
- accessories are fixed on the wooden top with height adjustment.

Accessories:

- horizontal spiral - for rotation in wrist joint,
- vertical roller with the spring - for dorsiflexion and palmar flexion in wrist joint,
- vertical ball with the spring - for abduction and adduction in wrist joint,
- pin - for improvement of manipulative skills,
- knob - for improvement of manipulative skills,
- arm stabilizer (stabilizer axis adjustment in the range of 0 - 90°).

Technical data:

Table top dimensions:	72 x 52 cm
Height adjustment:	54-86 cm
Mass:	8,4 kg