



Platform for gait and balance improvement PDB-2 with sandals

Link to product:

<http://meden.com.pl/offer/limb-dysfunction-therapy/80-platform-for-gait-and-balance-improvement-pdb-2-with-sandals.html>

Available online only

Condition: New

Informacje:

Hand crafted wobble balance platforms are quality exercise devices for balance and proprioceptive exercises. Designed to improve coordination, concentration and motor function, actively strengthen core muscles, ankles, knees, legs and back.

Product description:

Perfect to stimulate the inner ear to develop static and dynamic balance and gravitational orientation in conjunction with eyeball movement. It develops and improves movement smoothness. Made from high-quality plywood with hardwood ball. Quality lacquer finishing for easy maintenance and long lasting use. Optional version with sandals and straps attached to the pad.

Technical data:

Balance balls:	2
Platform shape:	rectangular
Platform size:	40 x 30 cm
Height:	72 mm
Max load:	150 kg