



Electrode tester Elza

Link to product:

<http://meden.com.pl/offer/other-medical-devices-and-equipment/19-electrode-tester-elza.html>

Available online only

Condition: New

Informacje:

Pocket-sized electrode tester for **electrotherapy**.

Product description:

Allows you for **quick** and **accurate examination** and cost reduction. The case has an **ergonomic shape** which allows you to conduct the test comfortably. The sound heard during the examination identifies suitability for further application of the tested electrode.

Interpretation of the results:

- **green LED indication** - efficient electrode
- **yellow-orange LED indicator** - partially drained electrode
- **red LED indicator** - totally drained or defective electrode

Technical data:

Power Supply: battery 12 V, typ 23 A