



Kriopol R - Bryza II

Link to product:


<http://meden.com.pl/offer/other-medical-devices-and-equipment/270-kriopol-r-bryza-ii.html>

Available online only

Condition: New

Informacje:

Cryotherapy device **Kriopol R - Bryza II**, is a product for cryotherapeutic procedures, often used by specialists in physiotherapy and cosmetic clinics.



Product features:

Category	Systemic cryotherapy
----------	----------------------

Product description:

Kriopol R - Bryza II

Benefits for patients:

- Reduced pain
- Increased metabolism
- Reduced activity of the inflammatory process
- Reduced use of anti-inflammatory and anaesthetic medications
- Reduced tension of skeletal muscles
- Faster regeneration and recovery
- Improved motion in the treated joints
- Faster treatment process
- Reduced swelling
- Alleviation of past-burn conditions

Advantages:

- Pressure-free, safe and easy to use
- Delivered with tank YDS-30 of capacity 31,5 l
- Different tank capacities YDS-15/30/50 corresponding to the client's needs (tank : 15 or 50 as option)
- Unique technology to maintain the elasticity of the tube supplying the nitrogen vapours onto the patient's skin
- Ready-to-use: 30 seconds after starting it reaches the working temperature of -160 degrees Celsius
- Does not cause cold burns, as the vapour stream is homogeneous and does not contain droplets of liquid nitrogen



Specification:

Temperature at nozzle exit [°C]:	-160
Temp. 10 cm, from nozzles [°C]:	-130
Time to reach working temperature:	30seconds
Display:	LED
Tank capacity [l]:	31,5 l (25,2 kg)
Number of treatments (YDS-30) per 1 filling (1 treatment = ca. 3 min.)	50-80
Dimensions (W x L x H)[cm]:	47 x 60 x 120
Blow adjustment steps:	4 + cryopuncture + 2 pulsed
Treatment time [min.]:	ca. 3
Coolant:	Liquid nitrogen
Power supply [V/Hz]:	230/50