



## Medical saddle stool with wheels FLEXI

Link to product:

<http://meden.com.pl/offer/other-medical-devices-and-equipment/32-medical-saddle-stool-with-wheels-flexi.html>

Available online only

Condition: New

### Informacje:

#### FLEXI Saddle stool:

**Lightweight** and **comfortable** stool for physicians and therapist with **chrome-plated** base and column



### Product description:

#### Features:

- poromeric leather upholstery
- chrome-plated base and column
- rubberised wheels
- angle-adjustment
- height adjustment
- upholstery available in colours matching the Prestige tables

#### Technical data:

<b>Height:</b>	500 - 700 mm
<b>Tilting adjustment</b>	15°
<b>Load weight:</b>	150 kg