



Urological and urodynamic examinations chair Mars

Link to product:

<http://meden.com.pl/offer/other-medical-devices-and-equipment/36-urological-and-urodynamic-examinations-chair-mars.html>

Available online only

Condition: New

Informacje:

A chair designed for gynaecology examination and treatment.



Product description:

Electric motors for easy adjustment of seat and backrest, including **Trendelenburg position**. Motors controlled with foot switch or **hand-operated remote controller**. Braked castors. Paper roll holder behind the backrest. Special bowl holder that folds under the seat. Optional panel turns the chair into a versatile operating table.

Main features:

- versatile
- ergonomic
- robust
- mobile
- smooth, precise adjustment
- high-quality finish

Standard accessories:

- leg rest (set of 2 pcs)
- arm rest (set of 2 pcs)
- basin with drainage
- funnel rim
- drip-bag hanger
- paper towel hanger
- manual and foot control
- rotary panels (R+L)
- table panel

Technical data:

Total length:	1820 mm
Total width:	900 mm
Width of seat:	600 mm
Width of back-rest:	600 mm

Minimum height of chair seat:	470 mm
Maximum height of chair seat:	1000 mm
Trendelenburg position:	20° ($\pm 2^\circ$)
Drive:	24 V
Max. capacity:	220 kg
Weight:	approx. 100
Power supply:	230 V ($\pm 10\%$) / 50 Hz