



Trolmed medical trolley type SPA-2

Link to product:

<http://meden.com.pl/offer/other-medical-devices-and-equipment/63-trolmed-medical-trolley-type-spa-2-.html>

Available online only

Condition: New

Informacje:

Trolmed trolley type SPA-2 designed mainly for physiotherapy equipment. Proper shelves configuration for convenient placement of various medical devices.



Product description:

Standard product:

- powder-painted aluminum frame
- 3 plastic shelves with raised edge preventing devices from slipping off
- four transport wheels (two front ones with brakes)

Technical data:

Depth:	43 cm
Width:	54 cm
Height:	89 cm
Number of shelves	1-3
Load capacity:	25 kg (shelf)