



Medical stool with wheels STANDARD

Upholstery colour Standard: blue,

Link to product:

<http://meden.com.pl/offer/other-medical-devices-and-equipment/64-medical-stool-with-wheels-standard.html>

Available online only

Condition: New

Informacje:

Round stool - Seat upholstered with poromeric leather.



Product description:

Standard:

- poromeric leather upholstery
- round seat
- base: chrome-plated steel
- castors: standard

Technical data:

Height:	520 - 720 mm
Seat diameter:	350 mm
Base diameter:	500 mm
Load weight:	150 kg