



Baby changing table SDN-1

Link to product:

<http://meden.com.pl/offer/other-medical-devices-and-equipment/69-baby-changing-table-sdn-1.html>

Available online only

Condition: New

Informacje:

SDN-1 baby changing table is your aid in care and examination of babies and infants at doctor's surgeries or at home.



Product description:

Optimal size work area allows for easy and safe baby care. Tall 5-inch (130 mm) sides for safety. Enables horizontal position in every condition by means of the adjustable table feet, guaranteeing the right level. Optionally, it can include a safe and handy side shelf and spacious bottom shelf for items such as scales, cosmetics, clothes, toys etc.

Robust, stable metal frame is covered with hygienic powder paint coating.

Padded surface upholstered with hygienic, easy-to-clean medical grade artificial leather. Upholstery has elegant, strong seams for long lasting use.