



1-section screen PDG-1

Link to product:

<http://meden.com.pl/offer/other-medical-devices-and-equipment/75-1-section-screen-pdg-1.html>

Available online only

Condition: New

Informacje:

Designed for surgeries, treatment rooms, patient rooms. Screen curtains made of medical grade fabrics (50 % cotton, 50% polyester) facilitating clean maintenance. **Powder - painted frame coating.** Steel construction guarantees **stability** and **durability**. In spite of the the steel framework the screen remains light.

Product description:

PDG-1 (1-section screen)

- steel framework
- two transport wheels

PDG-1

Width of the folded screen: 1000 mm

Width of the unfolded screen: 1000 mm

Height: 1805 mm