



Levitop Standard - Full Body Patient Floor Lift

Link to product:

<http://meden.com.pl/offer/patient-care/210-levitop-standard-full-body-patient-floor-lift.html>

Condition: New


Informacje:

Levitop Standard is full body patient floor lift are dedicated to patients who are temporarily or permanently unable to stand without assistance. Used to transfer them from bed to chair, wheelchair, chair as well as from the floor in case of falls.



Product description:

Levitop Standard - product features:





Hanger bar lock - in version Standard



Remote control - ergonomic and easy to use



Emergency stop button



Sling 4-point attachment - 4 hooks for loop sling fastening



2 castors with brake - brakes on the rear wheels prevent the device from moving during use when necessary

Technical data:

Lenght [cm]:	130
Width [cm]:	67
Height [cm]:	135
Base width (min.) [cm]:	52 / 67
Base width (max.) [cm]:	103 / 118
Lifting range [cm]:	68 / 196
Base legs height [cm]:	11,5
Weight [kg]:	61
Max. lifting capacity [kg]:	230
360° Rotating hanger:	V
Battery capacity [Ah]:	2,9
Charger [V/Hz]:	100-240 / 50 / 60 (integrated with control box)
Emergency stop button/ manual lowering mechanism:	V
Remote control functions:	<ul style="list-style-type: none"> • up / down • el. leg span adjustment • battery charge status <ul style="list-style-type: none"> • service indicator • overload indicator
Lifting speed without load (rising/lowering) [sec.]:	36 / 34
Turning radius (min./max.) [cm]:	144 / 158
Standard set:	<ul style="list-style-type: none"> • lift
Additional accessories:	<ul style="list-style-type: none"> • power cord for internal charger • manual sling to chose

Watch product presentation: