



Massage and treatment table Terapeuta Prestige M-P3

Link to product:

<http://meden.com.pl/offer/prestige-treatment-tables/85-massage-and-treatment-table-terapeuta-prestige-m-p3f0.html>

Available online only

Condition: New

Informacje:

3-section massage table.

Product description:

Product features:

- 3-section **massage and treatment table** with postural drainage setting
- Postural drainage with electric motor operated by 2 switches from both sides of the table
- Adjustable head section from +40° to -70° with gas spring
- Face hole on head section with removable face hole cover
- Special Aluminium Sliding System for augmented stability of leg section
- Adjustable leg section from 0° to 85° with gas spring
- Release handle of gas spring system in leg section accessible from both sides of the table
- Two-layer padding with non-flammable, bio-compatible and scratchproof upholstery in many colours
- Powder coated, extremely stable steel frame
- 2 side brackets for fixation of stabilization belts and elastic bands
- Traction support sleeves
- Personal Authorization System with 2 safety magnetic keys
- HalloTronic® System for electric high/low height adjustment from 51 to 99 cm with a capacity up to 225 kg
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table
- Non-slip rubber feet with 1 cm height adjustment allowing levelling of the table
- 2 types of base:
 - 4 castors with central lift system and 4 feet
 - 2 feet and 2 small wheels

Optional accessories:

- traction supports
- stabilization belt
- Kaltenborn wedge
- wedges
- bolsters
- holder for paper towel

Technical data:

Dimensions:	2050 x 690 mm
Number of sections:	3
Head-rest inclination angle:	+40/-70
Middle section inclination angle:	0/+30
Foot-rest inclination angle:	+85/0

Weight F0/F4:	95/110 kg
Maximum capacity:	225 kg
Power supply:	230 V / 50 Hz
Height adjustment:	51 - 99 cm