



## Massage and treatment table Terapeuta Prestige M-P7

Link to product:

<http://meden.com.pl/offer/prestige-treatment-tables/87-massage-and-treatment-table-terapeuta-prestige-m-p7.html>

Available online only

Condition: New

### Informacje:

7-section massage table.

### Product description:

#### Product features:

- 7 section massage and treatment table with postural drainage setting
- Adjustable head section from -70° to +40 with gas spring
- Face hole on head section with removable face hole cover
- Fold-down arm rests in head section
- Fold-down side sections in the middle section with semi-automatic release/blocking system
- Adjustable leg section from 0° to 85° with gas spring
- Special Aluminium Sliding System for augmented stability of leg section
- Release handle of gas spring system in leg section accessible from both sides of the table
- Two-layer padding with non-flammable, bio-compatible and scratchproof upholstery in many colours
- Powder coated, extremely stable steel frame
- 2 side brackets for fixation of stabilization belts and elastic bands
- Traction support sleeves
- Personal Authorization System with 2 safety magnetic keys
- HalloTronic® System for electric high/low height adjustment from 51 to 99 cm with a capacity up to 225 kg
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table
- Non-slip rubber feet with 1 cm height adjustment allowing levelling of the table
- 2 types of base:
  - basic wheeling system composed of 2 small wheels and 2 legs with non-slip rubber feet.
  - dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with nonslip rubber feet.

#### Optional accessories:

- traction supports
- stabilization belt
- Kaltenborn wedge
- wedges
- bolsters
- holder for paper towel

#### Technical data:

<b>Dimensions:</b>	2050 mm x 690 mm
<b>Number of sections:</b>	7
<b>Head-rest inclination angle:</b>	-70/+40
<b>Middle section inclination angle:</b>	0/+30

<b>Foot-rest inclination angle:</b>	0/+85
<b>Weight F0/F4:</b>	105/120 kg
<b>Maximum capacity:</b>	225 kg
<b>Power supply:</b>	230 V / 50 Hz
<b>Height adjustment:</b>	51 - 99 cm