



Massage and treatment table Safari Puma E-S2

Link to product:

<http://meden.com.pl/offer/safari-treatment-tables/189-massage-and-treatment-table-safari-puma-e-s2.html>

Available online only

Condition: New

Informacje:

2-section massage and treatment table with long head section

Product description:

Features:

- 2-section massage and treatment table with a long head section
- Width of table top: 69 or 80 cm
- Adjustable long head section from -20° to +85° or 0° to +85° (option) with 2 gas springs
- Round face hole in the head section with a removable face hole cover (2 types: round and crescent)
- Two-layer padding with non-flammable, bio-compatible and scratchproof upholstery in many colours
- Powder coated, extremely stable steel frame
- 2 side brackets for the fixation of stabilization belts and elastic bands
- Personal Authorization Safety System with 2 safety magnetic keys
- HalloTronic® System for electric high/low height adjustment from 52 to 93 cm with a capacity up to 200 kg
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table
- Non-slip rubber feet with 1 cm height adjustment allowing for the levelling of the table

Additional option:

- Upholstery of 80 mm thickness
- Special frame blocking the head section from going below 0°

Technical data:

Length [cm]:	211
Height [cm]:	52-93
Width [cm]:	69(80)
Weight F0/F4 [kg]	88/103
Max. load [kg]:	200

"I've been using several types of Meden-Inmed's Tables for over 15 years and I recommend them. I value their quality and design. I appreciate the care of the manufacturer about ergonomics of my work. I like the flexibility they give with all the adjustments they have. The lack of noise and stable work of motors and gas springs give my patients and me comfort. I recommend Meden-Inmed's Tables.

Lech Wojdyła
Neurac Instructor

Ordering information:

E-S2.F0: 2-section massage and treatment table with basic wheeling system composed of 2 small wheels and 2 legs with non-slip rubber feet.

E-S2.F4: 2-section massage and treatment table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with nonslip rubber feet.