

Nexus-H

Massage and treatment table



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**

Length: **202 cm**

Height: **48-105 cm**

Width: **70 cm**

Weight: **85 kg**

Max load: **185 kg**

Massage and treatment table **Nexus-H** is a perfect solution for manual therapy treating and mitigating injuries or impairments by restoring or improving massaged tissues and systems functions.

Massage and treatment table **Nexus-H** is based on a metal X-shaped base with excellent usability. It has been designed especially for masseurs and manual therapists, providing them stable and convenient table for their needs.

Using the latest technology, **Nexus-H** has matchless and eye-catching design and finish. Each of the components of the table is made with maximum precision, aesthetics and reliability.

With newest manufacturing technology and the most stringent production standards a stable, robust product was created, unrivalled in its price category.

Dimensions and operating parameters of the **Nexus-H** table have been developed with the help of engineers and specialists in the field of rehabilitation, cosmetology and massage. Thanks to the standards applied, we obtained a fully ergonomic design that provides comfort and safety for both the patient and therapist.

Features:

- Stable X-shaped base
- Durable construction with maximum load up to 185 kg
- Equipped with access authorization system to protect against unwanted table settings change
- Aesthetic and modern design perfectly fits the design of most rehabilitation and SPA centres
- Upholstery is made of durable, waterproof materials
- Wide range of upholstery colours
- Available in two and four section version

Applications:

- Rehabilitation clinics
- Beauty salons

Technical data:

Number of section:	2/4
Front section angle adjustment [°]:	-70°/ +40°
Height adjustment: [cm]:	48-105
Dimensions (length x width) [cm]:	202 x 70
Weight [kg]:	80
Maximum load [kg]:	185

