

# Matador

Massage and treatment table



Box dimensions: **215x 85 x 75 cm** Weight: **+25 kg**



Length: **202(205) cm**



Height: **48-95 cm**



Width: **69 cm**



Weight: **60 kg**



Max load: **160 kg**

Massage and treatment table **Matador** is a perfect choice for massage, kinesiotherapy and physical therapy.

Lightweight, durable, powder-painted steel frame with the upholstery made from top-quality artificial leather in many colours ensures reliability and convenience in therapist's daily work.

Massage and treatment table **Matador** is an excellent combination of high quality and low price. This makes it good solution for massage parlours, rehabilitation centers and physical therapy units.

## Features:

- Electric height adjustment with manual remote control
- Adjustable headrest from -70° to +40° with a removable face hole cover
- Non-slip rubber feet and 2 small wheels for moving of Matador table by lifting it from one side
- Non-flammable, bio-compatible and scratchproof upholstery of highest quality and functionality
- Wide range of upholstery colours
- Lightweight and stable construction
- Available in 1-, 2- or 3-section version

## Technical data:

|                                    |               |
|------------------------------------|---------------|
| Number of sections:                | 1/2/3         |
| Head section angle adjustment [°]: | -70°/+40°     |
| Height adjustment [cm]:            | 48-95         |
| Dimensions (length x width) [cm]:  | 202(205) x 69 |
| Weight [kg]:                       | 60            |
| Maximum load [kg]:                 | 160           |
| Paper roll width [cm]:             | 60            |
| Power supply [V/Hz]:               | 230/50        |

## Additional accessories:

- Attachable headrest (only 1-section version)
- Lifter for F4 frame
- Paper roll holder

